



ADVACRYL

Braided Coated Polyglactin 910 (Undyed & Violet)
Absorbable Surgical Suture U.S.P. (Synthetic)

**Consistent performance for confidence
in a variety of soft tissue approximation**



**ADVATRAY™ - the confidence to
suture with ease**



With great versatility comes the confidence to meet specific needs

Product Description

Sr. No.	Code	Suture			Needles		Units Per Box
		Colour	Size	Length (cm)	Description	Speciality	
1	ADVA2347WN180	VIOLET	1	180	-	-	12
2	ADVA2317WN180	VIOLET	2-0	180	-	-	12
3	ADVA2604MS	VIOLET	4-0	3 X 45	-	-	12

Sr. No.	Code	Suture			Needle		Units Per Box
		Colour	Size	Length (cm)	Description	Speciality	
4	ADVA2352L	VIOLET	2	110	40 mm, 1/2 Circle, Round Bodied	HEAVY	12
5	ADVA2452	VIOLET	2	90	40 mm, 1/2 Circle, Round Bodied	HEAVY	12
6	ADVA2353	VIOLET	2	90	50 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
7	ADVA2351	VIOLET	2	90	40 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
8	ADVA2478	VIOLET	2	90	40 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
9	ADVA2448	VIOLET	1	100	50 mm, 1/2 Circle, Round Bodied	HEAVY	12
10	ADVA2448L	VIOLET	1	110	50 mm, 1/2 Circle, Round Bodied	HEAVY	12
11	ADVA2326	VIOLET	1	70	45 mm, 3/8 Circle, Round Bodied	-	12
12	ADVA2647	VIOLET	1	90	45 mm, 1/2 Circle, Round Bodied	HEAVY	12
13	ADVA2347PB	VIOLET	1	100	45 mm, 1/2 Circle, Round Bodied	HEAVY	12
14	ADVA2344	UNDYED	1	90	40 mm, 1/2 Circle, Round Bodied	HEAVY	12
15	ADVA2347	VIOLET	1	90	40 mm, 1/2 Circle, Round Bodied	HEAVY	12
16	ADVA2347S	VIOLET	1	45	40 mm, 1/2 Circle, Round Bodied	HEAVY	12
17	ADVA2347VM	VIOLET	1	110	40 mm, 1/2 Circle, Round Bodied	HEAVY	12
18	ADVA2359	VIOLET	1	70	40 mm, 1/2 Circle, Round Bodied	HEAVY	12
19	ADVA2360	VIOLET	1	90	40 mm, 1/2 Circle, Round Bodied	-	12
20	ADVA2360L	VIOLET	1	120	40 mm, 1/2 Circle, Round Bodied	-	12
21	ADVA947	UNDYED	1	90	36 mm, 1/2 Circle, Round Bodied	-	12
22	ADVA2350	VIOLET	1	120	36 mm, 1/2 Circle, Round Bodied	-	12
23	ADVA2438	VIOLET	1	90	50 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
24	ADVA2421	VIOLET	1	90	40 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
25	ADVA2421G	VIOLET	1	90	40 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
26	ADVA2421S	VIOLET	1	45	40 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
27	ADVA2422	VIOLET	1	90	40 mm, 3/8 Circle, Reverse Cutting	-	12
28	ADVA2827	VIOLET	1	45	26 mm, 5/8 Circle, Reverse Cutting	-	12
29	ADVA2826	VIOLET	1	35	23 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
30	ADVA2826E	VIOLET	1	70	23 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
31	ADVA2364	VIOLET	1	90	40 mm, 1/2 Circle, Taper Cut	HEAVY	12
32	ADVA2363	VIOLET	1	70	40 mm, 1/2 Circle, Taper Cut	HEAVY	12

Sr. No.	Code	Suture			Needle		Units Per Box
		Colour	Size	Length (cm)	Description	Speciality	
33	ADVA2519	VIOLET	1	90	40 mm, 1/2 Circle, Taper Cut	-	12
34	ADVA2319	VIOLET	1	90	36 mm, 1/2 Circle, Taper Cut	-	12
35	ADVA9335	VIOLET	1	75	31 mm, J Type, Taper Cut	HEAVY	12
36	ADVA2511	VIOLET	1	45	30 mm, 1/2 Circle, Taper Cut	HEAVY	12
37	ADVA2346	VIOLET	0	90	40 mm, 1/2 Circle, Round Bodied	-	12
38	ADVA2346MM	VIOLET	0	110	40 mm, 1/2 Circle, Round Bodied	-	12
39	ADVA2346S	VIOLET	0	45	40 mm, 1/2 Circle, Round Bodied	-	12
40	ADVA2358	VIOLET	0	90	40 mm, 1/2 Circle, Round Bodied	-	12
41	ADVA2546	VIOLET	0	70	40 mm, 1/2 Circle, Round Bodied	-	12
42	ADVA946	UNDYED	0	90	36 mm, 1/2 Circle, Round Bodied	-	12
43	ADVA2337	UNDYED	0	90	30 mm, 1/2 Circle, Round Bodied	-	12
44	ADVA2338	VIOLET	0	90	30 mm, 1/2 Circle, Round Bodied	-	12
45	ADVA2342	VIOLET	0	70	30 mm, 1/2 Circle, Round Bodied	-	12
46	ADVA2533	VIOLET	0	90	40 mm, 3/8 Circle, Reverse Cutting	-	12
47	ADVA2534	VIOLET	0	90	36 mm, 1/2 Circle, Reverse Cutting	-	12
48	ADVA2511SB	VIOLET	0	45	26 mm, 5/8 Circle, Reverse Cutting	-	12
49	ADVA2825	VIOLET	0	45	23 mm, 1/2 Circle, Reverse Cutting	HEAVY	12
50	ADVA2510SB	VIOLET	0	45	23 mm, 1/2 Circle, Reverse Cutting	-	12
51	ADVA9901	UNDYED	0	75	40 mm, 1/2 Circle, Taper Cut	-	12
52	ADVA2518	VIOLET	0	90	40 mm, 1/2 Circle, Taper Cut	-	12
53	ADVA2518L	VIOLET	0	110	40 mm, 1/2 Circle, Taper Cut	-	12
54	ADVA1045	VIOLET	0	45	40 mm, 1/2 Circle, Taper Cut	-	12
55	ADVA2517	VIOLET	0	90	36 mm, 1/2 Circle, Taper Cut	-	12
56	ADVA2357	VIOLET	2-0	90	40 mm, 1/2 Circle, Round Bodied	-	12
57	ADVA2345	VIOLET	2-0	90	40 mm, 1/2 Circle, Round Bodied	-	12
58	ADVA945	UNDYED	2-0	90	36 mm, 1/2 Circle, Round Bodied	-	12
59	ADVA2317	VIOLET	2-0	90	30 mm, 1/2 Circle, Round Bodied	-	12
60	ADVA2317S	VIOLET	2-0	45	30 mm, 1/2 Circle, Round Bodied	-	12
61	ADVA2341	VIOLET	2-0	70	30 mm, 1/2 Circle, Round Bodied	-	12
62	ADVA2356	VIOLET	2-0	90	26 mm, 1/2 Circle, Round Bodied	HEAVY	12
63	ADVA602	VIOLET	2-0	70	26 mm, 5/8 Circle, Round Bodied	-	12
64	ADVA2356S	VIOLET	2-0	35	26 mm, 1/2 Circle, Round Bodied	-	12
65	ADVA2382	VIOLET	2-0	90	40 mm, 1/2 Circle, Reverse Cutting	-	12
66	ADVA2382S	VIOLET	2-0	45	40 mm, 1/2 Circle, Reverse Cutting	-	12
67	ADVA2316	VIOLET	2-0	90	30 mm, 3/8 Circle, Reverse Cutting	-	12
68	ADVA2694	VIOLET	2-0	45	26 mm, 5/8 Circle, Reverse Cutting	-	12
69	ADVA2404	UNDYED	2-0	90	26 mm, 3/8 Circle, Cutting	-	12

Sr. No.	Code	Suture			Needle		Units Per Box
		Colour	Size	Length (cm)	Description	Speciality	
70	ADVA945V	UNDYED	2-0	90	36 mm, 1/2 Circle, Taper Cut	-	12
71	ADVA2511A	VIOLET	2-0	45	30 mm, 1/2 Circle, Taper Cut	-	12
72	ADVA2511AA	VIOLET	2-0	20	30 mm, 1/2 Circle, Taper Cut	-	12
73	ADVA9350	VIOLET	2-0	75	26 mm, 1/2 Circle, Taper Cut	-	12
74	ADVA9828	VIOLET	2-0	75	22 mm, 1/2 Circle, Taper Cut	-	12
75	ADVA2122	VIOLET	2-0	70	26 mm, 1/2 Circle, Oval Round Bodied	-	12
76	ADVA322	VIOLET	3-0	70	36 mm, 1/2 Circle, Round Bodied	LONG TAPER	12
77	ADVA2528	VIOLET	3-0	90	25 mm, 1/2 Circle, Round Bodied	LONG TAPER	12
78	ADVA2437	VIOLET	3-0	70	20 mm, 1/2 Circle, Round Bodied	-	12
79	ADVA2427	VIOLET	3-0	70	12 mm, 1/2 Circle, Round Bodied	LONG TAPER	12
80	ADVA2407	UNDYED	3-0	70	26 mm, 3/8 Circle, Reverse Cutting	-	12
81	ADVA683	UNDYED	3-0	45	24 mm , 3/8 Circle, Reverse Cutting	PRIME	12
82	ADVA936	UNDYED	3-0	70	24 mm, 3/8 Circle, Reverse Cutting	-	12
83	ADVA2837	VIOLET	3-0	45	20 mm, 3/8 Circle, Reverse Cutting	-	12
84	ADVA2495	VIOLET	3-0	45	12 mm, 3/8 Circle, Reverse Cutting	PRIME	12
85	ADVA2328	UNDYED	3-0	90	25 mm, 3/8 Circle, Cutting	PRIME	12
86	ADVA2936	UNDYED	3-0	70	24 mm, 3/8 Circle, Cutting	-	12
87	ADVA2401	VIOLET	3-0	45	22 mm, 3/8 Circle, Cutting	-	12
88	ADVA2472	VIOLET	3-0	90	22 mm, 1/2 Circle, Cutting	-	12
89	ADVA2402	UNDYED	3-0	90	19 mm, 3/8 Circle, Cutting	PRIME	12
90	ADVA2402S	UNDYED	3-0	70	19 mm, 3/8 Circle, Cutting	PRIME	12
91	ADVA2471	VIOLET	3-0	45	16 mm, 3/8 Circle, Cutting	PRIME	12
92	ADVA2516	VIOLET	3-0	90	36 mm, 1/2 Circle, Taper Cut	-	12
93	ADVA2515	VIOLET	3-0	45	17 mm, 1/2 Circle, Taper cut	-	12
94	ADVA2123	VIOLET	3-0	70	26 mm, 1/2 Circle, Oval Round Bodied	-	12
95	ADVA2304	VIOLET	4-0	70	20 mm, 1/2 Circle, Round Bodied	-	12
96	ADVA2304S	VIOLET	4-0	45	20 mm, 1/2 Circle, Round Bodied	-	12
97	ADVA9106	VIOLET	4-0	75	17 mm, 1/2 Circle, Round Bodied	LONG TAPER	12
98	ADVA2305	VIOLET	4-0	70	16 mm, 1/2 Circle, Round Bodied	-	12
99	ADVA2494	VIOLET	4-0	45	12 mm, 3/8 Circle, Reverse Cutting	PRIME	12
100	ADVA2443	UNDYED	4-0	70	16 mm, 3/8 Circle, Cutting	PRIME	12
101	ADVA2443S	UNDYED	4-0	45	16 mm, 3/8 Circle, Cutting	PRIME	12
102	ADVA2465	VIOLET	4-0	45	16 mm, 5/8 Circle, Cutting	-	12
103	ADVA2403	VIOLET	5-0	45	16 mm, 1/2 Circle, Round Bodied	LONG TAPER	12
104	ADVA2395	VIOLET	5-0	70	10 mm, 3/8 Circle, Round Bodied	LONG TAPER	12
105	ADVA2872	UNDYED	5-0	45	13 mm, 3/8 Circle, Reverse Cutting	PRIME	12
106	ADVA2493	VIOLET	5-0	45	12 mm, 3/8 Circle, Reverse Cutting	PRIME	12
107	ADVA511G	UNDYED	5-0	45	10 mm, 1/2 Circle, Reverse Cutting	-	12
108	ADVA2442	VIOLET	5-0	45	16 mm, 3/8 Circle, Cutting	PRIME	12
109	ADVA2303	VIOLET	5-0	45	16 mm, 1/2 Circle, Oval Round Bodied	-	12

Sr. No.	Double Arm Code	Suture			Needle (1 & 2)		Units Per Box
		Colour	Size	Length (cm)	Description	Speciality	
110	ADVA2344DA180	VIOLET	1	180	40 mm, 1/2 Circle, Round Bodied (H) & 50 mm, 1/2 Circle, Round Bodied (H)	-	12
111	ADVA2347DA140	VIOLET	1	140	40 mm, 1/2 Circle, Round Bodied (H) & 40 mm, 1/2 Circle, Taper Cut	-	12
112	ADVA2347DA160	VIOLET	1	160	40 mm, 1/2 Circle, Round Bodied (H) & 40 mm, 1/2 Circle, Taper Cut	-	12
113	ADVA2347DA180	VIOLET	1	180	40 mm, 1/2 Circle, Round Bodied (H) & 40 mm, 1/2 Circle, Taper Cut	-	12
114	ADVA2347DA180H	VIOLET	1	180	40 mm, 1/2 Circle, Round Bodied (H) & 40 mm, 1/2 Circle, Taper Cut	-	36
115	ADVA2346DA180	VIOLET	0	180	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut	-	12
116	ADVA2318	VIOLET	2-0	90	30 mm, 1/2 Circle, Round Bodied	-	12
117	ADVA2436	VIOLET	3-0	90	26 mm, 1/2 Circle, Round Bodied	LONG TAPER	12
118	ADVA2354	VIOLET	4-0	90	16 mm, 3/8 Circle, Cutting PRIME & 16 mm, 1/2 Circle, Round Bodied	-	12
119	ADVA2670	VIOLET	6-0	45	8 mm, 1/4 Circle, Spatulated Micro Point	-	12

Product Description


Composition	90% Glycolide, 10% L-lactide
Coating	Copolymer of Poly (glycolide-co-lactide) and Calcium Stearate
Size	U.S.P. size 2 to 6-0
Colour	Undyed and Violet
Absorption Profile	Essentially complete up to 80 days
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Trocar, Reverse Cutting, Taper Cut, Blunt Point, Oval Round Bodied, Diamond Point & Spatula Micro Point needles

Tensile Strength

Days of Implantation	Tensile Strength Retention
14 Days	70%
21 Days	50%
28 Days	20%

34% more knot pull tensile strength than U.S.P. limit ensures knot security


Knot pull tensile strength (KgF)¹



REFERENCE _____

1 : Data on file with AMS - Bench test done on Polyglactin-910 U.S.P. size - 1, 40 mm, Round bodied needle


Multifilament Material



Provides excellent support for a range of critical wound healing period to secure closure in a variety of tissues, procedures and patients.

ADVA TRAY™ - the confidence to suture with ease

ADVA TRAY™ a unique dispensing system designed to deliver ease of suture dispensing with minimal memory in the suture strand



Premium Needles for consistent Superior Performance

PREMIUM NEEDLES
ensure
UNSURPASSED SUTURING EXPERIENCE

DURABLE COAT RETAINS SHARPNESS... PASS AFTER PASS	Maintains needle sharpness over the duration of procedure, resulting into consistent penetration pass after pass thus less potential for tissue trauma.
RIBBED BODY FOR SUPERIOR CONTROL	For superior grasping and stability of needle in needle forceps.
CHAMFERING MINIMUM SUTURE DETACHMENT	Diminishes chances for future fracturing at the swage, by rounding off the sharp edges of the barrel in the swage area of the needle.
MICRO-FINISHING SMOOTH & BRIGHT	Makes needle surface free from micro asperities, thus creates ultra sharp edges for lesser tissue drag.

Disclaimer – These features are not common across the needle families

Minimal tissue trauma enables faster recovery



Advanced MedTech Solutions

Advanced MedTech Solutions Private Limited

Corporate : 305, Ceejay House, F-Block, Shiv Sagar Estate, Dr. Annie Besant Road, Worli, Mumbai 400018, India.

T : +91 22 4946100 / 01 | F : +91 22 24946102 | E: info@amsltd.com | amsltd.com

Plant : P-21-22, 25-26 & 34-35, GIDC Manjusar, Tal - Savli, Dist. - Vadodara - 391775 Gujarat, India | T: +91 2667 671302

Advanced MedTech Solutions Pte. Ltd.

111 Somerset Road, #07-08 111 Somerset, Singapore - 238164 | T: (65) 636910 80