



ADVASYL

Proteinaceous silk fibres called fibroin
Braided Natural Non-Absorbable Surgical Suture

Exceptional control for confidence in
gentle tissue approximation





Excellent handling combined with great strength for gentle tissue approximation

Product Description

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVAW075	BLACK	0	75	-	-	-	-	12
2	ADVA204	BLACK	0	17 X 45	-	-	-	-	12
3	ADVA214	BLACK	0	2 X 76	-	-	-	-	12
4	ADVA224	BLACK	0	10 X 35	-	-	-	-	12
5	ADVA234	BLACK	0	2 X 100	-	-	-	-	12
6	ADVA5109	BLACK	0	13 X 60	-	-	-	-	12
7	ADVA215	BLACK	1	2 X 76	-	-	-	-	12
8	ADVA225	BLACK	1	10 X 75	-	-	-	-	12
9	ADVA5108	BLACK	1	180	-	-	-	-	12
10	ADVA5121	BLACK	1	17 X 45	-	-	-	-	12
11	ADVA5155	BLACK	1	13 X 60	-	-	-	-	12
12	ADVA216	BLACK	2	2 X 76	-	-	-	-	12
13	ADVA216L	BLACK	2	180	-	-	-	-	12
14	ADVA5165	BLACK	2	12 X 75	-	-	-	-	12
15	ADVA213	BLACK	2-0	2 X 76	-	-	-	-	12
16	ADVAW218	BLACK	2-0	180	-	-	-	-	12
17	ADVA223	BLACK	2-0	10 X 35	-	-	-	-	12
18	ADVA5122	BLACK	2-0	17 X 45	-	-	-	-	12
19	ADVA5136	BLACK	2-0	13 X 60	-	-	-	-	12
20	ADVA5138	BLACK	2-0	15 X 60	-	-	-	-	12
21	ADVA5169	BLACK	2-0	75	-	-	-	-	12
22	ADVA212	BLACK	3-0	2 X 76	-	-	-	-	12
23	ADVA222	BLACK	3-0	10 X 35	-	-	-	-	12
24	ADVA232	BLACK	3-0	2 X 100	-	-	-	-	12
25	ADVAW318	BLACK	3-0	180	-	-	-	-	12
26	ADVA5106	BLACK	3-0	13 X 60	-	-	-	-	12
27	ADVA5107	BLACK	3-0	17 X 45	-	-	-	-	12
28	ADVA5140	BLACK	3-0	15 X 60	-	-	-	-	12
29	ADVA5168	BLACK	3-0	60	-	-	-	-	12
30	ADVA5187	BLACK	3-0	30 X 75	-	-	-	-	12

⊙ Round Bodied

▽ Reverse Cutting

△ Cutting

⊕ Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
31	ADVA211	BLACK	4-0	2 X 76	-	-	-	-	12
32	ADVAW418	BLACK	4-0	180	-	-	-	-	12
33	ADVA5123	BLACK	4-0	13 X 60	-	-	-	-	12
34	ADVA5167	BLACK	4-0	60	-	-	-	-	12
35	ADVA5186	BLACK	4-0	30 X 75	-	-	-	-	12
36	ADVA334	BLACK	0	75	1/2 CIRCLE		31	-	12
37	ADVA334H	BLACK	0	75	1/2 CIRCLE		31	-	36
38	ADVA782H	BLACK	0	75	3/8 CIRCLE		45	-	36
39	ADVA5016	BLACK	0	75	3/8 CIRCLE		36	-	12
40	ADVA5017	BLACK	0	76	3/8 CIRCLE		30	-	12
41	ADVA5037	BLACK	0	76	3/8 CIRCLE		45	-	12
42	ADVA5037L	BLACK	0	90	3/8 CIRCLE		45	-	12
43	ADVA5038	BLACK	0	76	1/2 CIRCLE		45	-	12
44	ADVA5110	BLACK	0	75	STRAIGHT		60	-	12
45	ADVA5130	BLACK	0	70	1/2 CIRCLE		40	-	12
46	ADVA5132	BLACK	0	75	3/8 CIRCLE		48	-	12
47	ADVA5134	BLACK	0	75	1/2 CIRCLE		40	-	12
48	ADVA5142	BLACK	0	75	1/2 CIRCLE		26	-	12
49	ADVA5147	BLACK	0	100	1/2 CIRCLE		90	-	12
50	ADVA5152	BLACK	0	75	3/8 CIRCLE		30	-	12
51	ADVA5156	BLACK	0	75	3/8 CIRCLE		40	-	12
52	ADVA5175	BLACK	0	90	1/2 CIRCLE		30	-	12
53	ADVA5176	BLACK	0	70	1/2 CIRCLE		30	-	12
54	ADVA5179	BLACK	0	70	1/2 CIRCLE		25	-	12
55	ADVA5180	BLACK	0	70	3/8 CIRCLE		45	-	12
56	ADVA5197	BLACK	0	75	1/2 CIRCLE		40	-	12
57	ADVA5203	BLACK	0	75	1/2 CIRCLE		36	-	12
58	ADVA5209	BLACK	0	75	3/8 CIRCLE		45	-	12
59	ADVA5279	BLACK	0	70	1/2 CIRCLE		25	-	12
60	ADVA5332	BLACK	0	90	1/2 CIRCLE		30	-	12
61	ADVA5334	BLACK	0	76	1/2 CIRCLE		30	-	12
62	ADVA5575	BLACK	0	75	3/8 CIRCLE		55	-	12
63	ADVA5750	BLACK	0	75	1/2 CIRCLE		50	-	12
64	ADVA8351545	BLACK	0	75	3/8 CIRCLE		45	-	12
65	ADVA723H	BLACK	1	75	1/2 CIRCLE		31	HEAVY	36
66	ADVA5005	BLACK	1	76	1/2 CIRCLE		38	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
67	ADVA5062	BLACK	1	76	3/8 CIRCLE		60	-	12
68	ADVA5062L	BLACK	1	90	3/8 CIRCLE		60	-	12
69	ADVA5062H	BLACK	1	76	3/8 CIRCLE		60	-	36
70	ADVA5063	BLACK	1	76	3/8 CIRCLE		50	-	12
71	ADVA5101	BLACK	1	75	1/2 CIRCLE		36	-	12
72	ADVA5101H	BLACK	1	75	1/2 CIRCLE		36	-	36
73	ADVA5102	BLACK	1	75	1/2 CIRCLE		36	-	12
74	ADVA5124	BLACK	1	75	1/2 CIRCLE		35	-	12
75	ADVA5129	BLACK	1	70	1/2 CIRCLE		40	-	12
76	ADVA5133	BLACK	1	75	3/8 CIRCLE		48	-	12
77	ADVA5150	BLACK	1	75	1/2 CIRCLE		60	-	12
78	ADVA5157	BLACK	1	75	3/8 CIRCLE		40	-	12
79	ADVA5158	BLACK	1	75	1/2 CIRCLE		45	-	12
80	ADVA5162	BLACK	1	90	1/2 CIRCLE		40	-	12
81	ADVA5171	BLACK	1	75	3/8 CIRCLE		24	-	12
82	ADVA5173	BLACK	1	76	3/8 CIRCLE		36	-	12
83	ADVA5177	BLACK	1	75	1/2 CIRCLE		30	-	12
84	ADVA5202	BLACK	1	75	3/8 CIRCLE		60	-	12
85	ADVA5634	BLACK	1	75	1/2 CIRCLE		30	-	12
86	ADVA8401545	BLACK	1	75	3/8 CIRCLE		45	-	12
87	ADVA5160	BLACK	2	90	1/2 CIRCLE		45	-	12
88	ADVA5161	BLACK	2	90	1/2 CIRCLE		45	-	12
89	ADVA5172	BLACK	2	75	1/2 CIRCLE		40	HEAVY	12
90	ADVA5178	BLACK	2	75	1/2 CIRCLE		30	-	12
91	ADVA5198	BLACK	2	75	1/2 CIRCLE		4	HEAVY	12
92	ADVA5061	BLACK	4	76	1/2 CIRCLE		35	-	12
93	ADVA5061SP	BLACK	4	76	1/2 CIRCLE		40	-	12
94	ADVA321	BLACK	2-0	45	3/8 CIRCLE		26	-	12
95	ADVA321H	BLACK	2-0	45	3/8 CIRCLE		26	-	36
96	ADVA333	BLACK	2-0	75	1/2 CIRCLE		31	-	12
97	ADVA333H	BLACK	2-0	75	1/2 CIRCLE		31	-	36
98	ADVA791	BLACK	2-0	100	3/8 CIRCLE		90	-	12
99	ADVA5036	BLACK	2-0	76	3/8 CIRCLE		45	-	12
100	ADVA5036L	BLACK	2-0	90	3/8 CIRCLE		45	-	12
101	ADVA5052	BLACK	2-0	76	1/2 CIRCLE		50	-	12
102	ADVA5065	BLACK	2-0	76	1/2 CIRCLE		26	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
103	ADVA5104	BLACK	2-0	76	1/2 CIRCLE		26	-	12
104	ADVA5111	BLACK	2-0	75	1/2 CIRCLE		35	-	12
105	ADVA5112	BLACK	2-0	100	STRAIGHT		55	-	12
106	ADVA5112H	BLACK	2-0	100	STRAIGHT		55	-	36
107	ADVA5125	BLACK	2-0	75	1/2 CIRCLE		25	-	12
108	ADVA5126	BLACK	2-0	75	1/2 CIRCLE		22	-	12
109	ADVA5127	BLACK	2-0	45	1/2 CIRCLE		26	-	12
110	ADVA5137	BLACK	2-0	75	3/8 CIRCLE		26	-	12
111	ADVA5139	BLACK	2-0	75	3/8 CIRCLE		35	-	12
112	ADVA5146	BLACK	2-0	45	1/2 CIRCLE		26	PRIME	12
113	ADVA5149	BLACK	2-0	75	1/2 CIRCLE		40	-	12
114	ADVA5151	BLACK	2-0	100	1/2 CIRCLE		90	-	12
115	ADVA5153	BLACK	2-0	75	3/8 CIRCLE		30	-	12
116	ADVA5159	BLACK	2-0	75	3/8 CIRCLE		24	-	12
117	ADVA5164	BLACK	2-0	90	3/8 CIRCLE		24	-	12
118	ADVA5170	BLACK	2-0	75	1/2 CIRCLE		25	-	12
119	ADVA5181	BLACK	2-0	70	3/8 CIRCLE		45	-	12
120	ADVA5182	BLACK	2-0	70	3/8 CIRCLE		45	-	12
121	ADVA5190	BLACK	2-0	75	STRAIGHT		60	-	12
122	ADVA5194	BLACK	2-0	45	1/2 CIRCLE		24	-	12
123	ADVA5199	BLACK	2-0	45	1/2 CIRCLE		26	-	12
124	ADVA5204	BLACK	2-0	45	3/8 CIRCLE		26	PRIME	12
125	ADVA5208	BLACK	2-0	75	3/8 CIRCLE		26	-	12
126	ADVA5210	BLACK	2-0	100	1/2 CIRCLE		40	-	12
127	ADVA5224	BLACK	2-0	75	1/2 CIRCLE		24	-	12
128	ADVA5236	BLACK	2-0	75	3/8 CIRCLE		36	-	12
129	ADVA5275	BLACK	2-0	75	1/2 CIRCLE		22	-	12
130	ADVA5331	BLACK	2-0	90	1/2 CIRCLE		30	-	12
131	ADVA5333	BLACK	2-0	76	1/2 CIRCLE		30	-	12
132	ADVA5611	BLACK	2-0	90	1/2 CIRCLE		17	-	12
133	ADVA5670	BLACK	2-0	76	1/2 CIRCLE		25	-	12
134	ADVA576	BLACK	3-0	75	5/8 CIRCLE		25	-	12
135	ADVA5205	BLACK	3-0	45	1/2 CIRCLE		22	PRIME	12
136	ADVA5206	BLACK	3-0	45	3/8 CIRCLE		26	PRIME	12
137	ADVA5002	BLACK	3-0	76	3/8 CIRCLE		16	-	12
138	ADVA5003	BLACK	3-0	90	3/8 CIRCLE		16	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
139	ADVA5028	BLACK	3-0	76	3/8 CIRCLE		26	-	12
140	ADVA5028S	BLACK	3-0	45	3/8 CIRCLE		26	-	12
141	ADVA5028SH	BLACK	3-0	45	3/8 CIRCLE		26	-	36
142	ADVA5064	BLACK	3-0	76	1/2 CIRCLE		22	-	12
143	ADVA5070	BLACK	3-0	76	1/2 CIRCLE		25	-	12
144	ADVA5085	BLACK	3-0	90	1/2 CIRCLE		25	-	12
145	ADVA5087	BLACK	3-0	76	1/2 CIRCLE		20	-	12
146	ADVA5100	BLACK	3-0	90	1/2 CIRCLE		35	-	12
147	ADVA5103	BLACK	3-0	76	3/8 CIRCLE		26	-	12
148	ADVA5113	BLACK	3-0	75	1/2 CIRCLE		16	-	12
149	ADVA5114	BLACK	3-0	75	1/2 CIRCLE		35	-	12
150	ADVA5120	BLACK	3-0	75	1/2 CIRCLE		16	-	12
151	ADVA5128	BLACK	3-0	45	1/2 CIRCLE		26	-	12
152	ADVA5131	BLACK	3-0	45	1/2 CIRCLE		26	-	12
153	ADVA5135	BLACK	3-0	75	1/2 CIRCLE		30	-	12
154	ADVA5141	BLACK	3-0	75	3/8 CIRCLE		19	-	12
155	ADVA5143	BLACK	3-0	75	3/8 CIRCLE		24	-	12
156	ADVA5144	BLACK	3-0	75	3/8 CIRCLE		17	-	12
157	ADVA5145	BLACK	3-0	45	1/2 CIRCLE		25	SLIM BLADE	12
158	ADVA5154	BLACK	3-0	75	3/8 CIRCLE		25	-	12
159	ADVA5166	BLACK	3-0	75	1/2 CIRCLE		40	-	12
160	ADVA5183	BLACK	3-0	70	3/8 CIRCLE		45	-	12
161	ADVA5184	BLACK	3-0	60	3/8 CIRCLE		22	-	12
162	ADVA5188	BLACK	3-0	45	3/8 CIRCLE		24	-	12
163	ADVA5189	BLACK	3-0	75	1/2 CIRCLE		25	-	12
164	ADVA5191	BLACK	3-0	45	3/8 CIRCLE		19	-	12
165	ADVA5193	BLACK	3-0	75	3/8 CIRCLE		24	-	12
166	ADVA5263	BLACK	3-0	75	1/2 CIRCLE		26	-	12
167	ADVA5263H	BLACK	3-0	75	1/2 CIRCLE		26	-	36
168	ADVA5319	BLACK	3-0	75	3/8 CIRCLE		19	-	12
169	ADVA5326	BLACK	3-0	75	1/2 CIRCLE		26	-	12
170	ADVA5326XL	BLACK	3-0	90	1/2 CIRCLE		26	-	12
171	ADVA6571H	BLACK	3-0	75	1/2 CIRCLE		31	-	36
172	ADVA329	BLACK	4-0	45	3/8 CIRCLE		19	-	12
173	ADVA329H	BLACK	4-0	45	3/8 CIRCLE		19	-	36
174	ADVA582	BLACK	4-0	75	3/8 CIRCLE		17	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
175	ADVA2505	BLACK	4-0	75	1/2 CIRCLE		20	-	12
176	ADVA5000	BLACK	4-0	90	3/8 CIRCLE		16	PRIME	12
177	ADVA5001	BLACK	4-0	76	3/8 CIRCLE		16	-	12
178	ADVA5049	BLACK	4-0	76	1/2 CIRCLE		16	-	12
179	ADVA5050	BLACK	4-0	90	1/2 CIRCLE		16	-	12
180	ADVA5082	BLACK	4-0	76	3/8 CIRCLE		16	LONG TAPER	12
181	ADVA5086	BLACK	4-0	76	1/2 CIRCLE		20	-	12
182	ADVA5086H	BLACK	4-0	76	1/2 CIRCLE		20	-	36
183	ADVA5099	BLACK	4-0	76	1/2 CIRCLE		17	-	12
184	ADVA5115	BLACK	4-0	45	1/2 CIRCLE		19	-	12
185	ADVA5116	BLACK	4-0	75	1/2 CIRCLE		16	-	12
186	ADVA5148	BLACK	4-0	45	1/2 CIRCLE		25	-	12
187	ADVA5185	BLACK	4-0	75	1/2 CIRCLE		25	-	12
188	ADVA5185XL	BLACK	4-0	90	1/2 CIRCLE		26	-	12
189	ADVA5192	BLACK	4-0	45	1/2 CIRCLE		16	-	12
190	ADVA5196	BLACK	4-0	75	3/8 CIRCLE		19	-	12
191	ADVA5201	BLACK	4-0	45	3/8 CIRCLE		25	SLIM BLADE	12
192	ADVA5207	BLACK	4-0	4 X 45	3/8 CIRCLE		16	-	12
193	ADVA5416	BLACK	4-0	75	3/8 CIRCLE		16	-	12
194	ADVA5416H	BLACK	4-0	75	3/8 CIRCLE		16	-	36
195	ADVA64086	BLACK	4-0	75	1/2 CIRCLE		13	-	12
196	ADVA500H	BLACK	5-0	75	3/8 CIRCLE		16	-	36
197	ADVA581	BLACK	5-0	75	3/8 CIRCLE		17	-	12
198	ADVA5027	BLACK	5-0	76	3/8 CIRCLE		12	LONG TAPER	12
199	ADVA5079	BLACK	5-0	90	3/8 CIRCLE		12	-	12
200	ADVA5080	BLACK	5-0	76	3/8 CIRCLE		12	-	12
201	ADVA5081	BLACK	5-0	76	3/8 CIRCLE		16	LONG TAPER	12
202	ADVA5081SP	BLACK	5-0	76	3/8 CIRCLE		16	PRIME	12
203	ADVA5095	BLACK	5-0	76	3/8 CIRCLE		10	LONG TAPER	12
204	ADVA5117	BLACK	5-0	75	1/2 CIRCLE		16	-	12
205	ADVA5195	BLACK	5-0	75	1/2 CIRCLE		19	-	12
206	ADVA5200	BLACK	5-0	45	3/8 CIRCLE		15	SLIM BLADE	12
207	ADVA5211	BLACK	5-0	75	3/8 CIRCLE		19	-	12
208	ADVA5315	BLACK	5-0	45	3/8 CIRCLE		13	-	12
209	ADVA5695	BLACK	5-0	76	3/8 CIRCLE		12	-	12
210	ADVA529	BLACK	6-0	45	3/8 CIRCLE		16	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
211	ADVA529H	BLACK	6-0	45	3/8 CIRCLE		16	-	36
212	ADVA681SP	BLACK	6-0	76	1/2 CIRCLE		16	-	12
213	ADVA5012	BLACK	6-0	38	3/8 CIRCLE		10	MICRO POINT	12
214	ADVA5029	BLACK	6-0	76	3/8 CIRCLE		13	-	12
215	ADVA5029SP	BLACK	6-0	76	3/8 CIRCLE		11	PRIME	12
216	ADVA5043	BLACK	6-0	38	3/8 CIRCLE		8	MICRO POINT	12
217	ADVA5088	BLACK	6-0	38	1/4 CIRCLE		8	MICRO POINT	12
218	ADVA5119	BLACK	6-0	45	1/2 CIRCLE		15	SLIM BLADE	12
219	ADVA5612	BLACK	6-0	45	3/8 CIRCLE		12	-	12
220	ADVA5019	BLACK	8-0	38	1/2 CIRCLE		8	MICRO POINT	12
221	ADVA5022	BLACK	8-0	38	3/8 CIRCLE		6	-	12
222	ADVA5174	BLACK	8-0	60	1/2 CIRCLE		13	-	12

Sr. No.	Double Arm Code	Suture			Needle (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
223	ADVA5118	BLACK	5-0	45	1/2 CIRCLE			10	12
224	ADVA888	BLACK	6-0	45	1/4 CIRCLE			8	12
225	ADVA1732	BLACK	6-0	45	3/8 CIRCLE			6.5	12
226	ADVA5105	BLACK	6-0	75	3/8 CIRCLE	MICRO POINT 	MICRO POINT 	12	12
227	ADVA5021	BLACK	8-0	38	3/8 CIRCLE			6	12

 Round Bodied

 Reverse Cutting

 Cutting

 Spatulated Micropoint

Product Description

Composition	Proteinaceous silk fibres called fibroin
Coating	Wax
Sizes	U.S.P. Size 4 to 6-0 & 8-0
Colour	Black
Absorption Profile	Non-Absorbable
Needle Availability	Available as Pre-cut sutures and in wide variety of Ball Point, Cutting, Round Bodied (Taper Point), Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied, Spatulated Micro Point & Diamond Point

Multifilament Material




Black

Tensile Strength

Permanent however gradual loss in tensile strength over a period of time.

60% more knot pull tensile strength Vs. U.S.P. limit ensures high tensile strength offering strong & secure closure for permanent wound support.


Knot pull tensile strength (KgF)¹



REFERENCE _____

1 : Data on file with AMS - Bench test done on Silk U.S.P. Size - 2-0, 45 mm Reverse cutting needle

Premium Needles for consistent Superior Performance



Minimal tissue trauma enables faster recovery



P-21-22, 25-26 & 34-35, GIDC, Manjusar, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India
T : (91) 2667 671302 | amsltd.com | info@amsltd.com

Singapore Office : 111, Somerset Road, 07-08, 111 Somerset, Singapore - 238164.
T: (65) 636910 80 / 89 | E: info@amsltd.com | amsltd.com