



ADVABOND

Polyester Green & White (Undyed)

Braided Coated Non Absorbable Surgical Suture

Ultimate security for confidence in approximation of tissues under stress





Enhancing confidence for everlasting support

Product Description

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVA6875	GREEN	0	75	1/2 CIRCLE		30	-	12
2	ADVA6336SN	GREEN	0	75	1/2 CIRCLE		36	-	12
3	ADVA26532	GREEN	0	75	3/8 CIRCLE		30	-	12
4	ADVA0969	GREEN	1	75	J Type		31	-	12
5	ADVA0538	GREEN	1	75	3/8 CIRCLE		36	-	12
6	ADVA6016	GREEN	1	90	1/2 CIRCLE		26	-	12
7	ADVA6445	GREEN	1	45	1/2 CIRCLE		45	-	12
8	ADVA643	GREEN	2	100	1/2 CIRCLE		45	HEAVY	12
9	ADVA4842	GREEN	2	75	1/2 CIRCLE		45	HEAVY	12
10	ADVA4843	GREEN	2	4 X 75	1/2 CIRCLE		45	HEAVY	12
11	ADVA6023	GREEN	2	90	1/2 CIRCLE		55	-	12
12	ADVA6014	GREEN	4	90	1/2 CIRCLE		50	-	12
13	ADVA645	GREEN	5	100	1/2 CIRCLE		55	HEAVY	12
14	ADVA0B46	GREEN	5	4 X 75	1/2 CIRCLE		48	-	12
15	ADVA4845	GREEN	5	75	1/2 CIRCLE		55	HEAVY	12
16	ADVA4846	GREEN	5	4 X 75	1/2 CIRCLE		55	HEAVY	12
17	ADVA4846L	GREEN	5	90	1/2 CIRCLE		55	-	12
18	ADVA6001	GREEN	5	100	3/8 CIRCLE		50	-	12
19	ADVA6015	GREEN	5	90	1/2 CIRCLE		50	-	6
20	ADVA26191	GREEN	5	75	1/2 CIRCLE		37	-	12
21	ADVA673A	WHITE	2-0	76	1/2 CIRCLE		17	-	12
22	ADVAX873	GREEN	2-0	75	1/2 CIRCLE		17	LONG TAPER	12
23	ADVA6002	GREEN	2-0	75	1/2 CIRCLE		26	-	12
24	ADVA6011	GREEN	2-0	100	1/2 CIRCLE		25	-	12
25	ADVA6887	GREEN	2-0	90	3/8 CIRCLE		25	-	6
26	ADVA6916	GREEN	2-0	90	1/2 CIRCLE		17	-	6
27	ADVA6917SP	WHITE	2-0	75	1/2 CIRCLE		17	-	12
28	ADVA26531	GREEN	2-0	75	3/8 CIRCLE		30	-	12
29	ADVA655	GREEN	3-0	100	1/2 CIRCLE		25	LONG TAPER	12
30	ADVA6007	GREEN	3-0	75	1/2 CIRCLE		26	-	12
31	ADVA6542	GREEN	3-0	75	1/2 CIRCLE		25	-	12

 Round Bodied

 Cutting

 Taper Cut

H - Heavy

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
32	ADVA6832	GREEN	3-0	75	1/2 CIRCLE		26	LONG TAPER	12
33	ADVA6836	GREEN	3-0	75	1/2 CIRCLE		17	-	12
34	ADVA6004	GREEN	4-0	45	1/2 CIRCLE		10	-	12
35	ADVA6513	GREEN	5-0	75	1/2 CIRCLE		13	-	12
36	ADVA26551	GREEN	5-0	45	3/8 CIRCLE		16	-	12
37	ADVA6010	GREEN	6-0	75	1/2 CIRCLE		13	-	12
38	ADVA6024	GREEN	2	180	-	-	-	-	12

Sr. No.	Double Arm Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
39	ADVA6008	GREEN	0	90	1/2 CIRCLE			26	12
40	ADVA6028	GREEN	0	90	1/2 CIRCLE			25	12
41	ADVA6336	GREEN	0	75	1/2 CIRCLE			36	12
42	ADVA4846DA	GREEN	5	75	1/2 CIRCLE	 H	 H	55	12
43	ADVA617	GREEN	2-0	90	1/2 CIRCLE			17	12
44	ADVA617A	WHITE	2-0	90	1/2 CIRCLE			17	12
45	ADVAP617	GREEN	2-0	90	1/2 CIRCLE			17	12
46	ADVAP617A	WHITE	2-0	90	1/2 CIRCLE			17	12
47	ADVAP617S	GREEN	2-0	75	1/2 CIRCLE			17	12
48	ADVA687	GREEN	2-0	90	1/2 CIRCLE			25	12
49	ADVA687A	WHITE	2-0	90	1/2 CIRCLE			25	12
50	ADVAP687	GREEN	2-0	90	1/2 CIRCLE			25	12
51	ADVAP687A	WHITE	2-0	90	1/2 CIRCLE			25	12
52	ADVAP687S	GREEN	2-0	75	1/2 CIRCLE			25	12
53	ADVAP687PL3	GREEN	2-0	90	1/2 CIRCLE			25	12
54	ADVAP687APL3	WHITE	2-0	90	1/2 CIRCLE			25	12
55	ADVA697	GREEN	2-0	75	1/2 CIRCLE			25	12
56	ADVA697A	WHITE	2-0	75	1/2 CIRCLE			25	12
57	ADVA6009	WHITE	2-0	90	1/2 CIRCLE			17	12
58	ADVA6012	GREEN	2-0	75	1/2 CIRCLE			18	12
59	ADVA6013PL3	GREEN	2-0	75	1/2 CIRCLE			20	12
60	ADVA6017PL6	GREEN	2-0	75	1/2 CIRCLE	 H	 H	20	12
61	ADVA6017APL6	WHITE	2-0	75	1/2 CIRCLE	 H	 H	20	12
62	ADVA6018	GREEN	2-0	90	1/2 CIRCLE			20	12
63	ADVA6018A	WHITE	2-0	90	1/2 CIRCLE			20	12
64	ADVA6019PL6	GREEN	2-0	90	1/2 CIRCLE			25	12
65	ADVA6019APL6	WHITE	2-0	90	1/2 CIRCLE			25	12
66	ADVA6025	GREEN	2-0	90	1/2 CIRCLE			17	12
67	ADVA6219	GREEN	2-0	90	1/2 CIRCLE			17	12
68	ADVA6225	GREEN	2-0	90	1/2 CIRCLE			25	12
69	ADVA6425	WHITE	2-0	90	1/2 CIRCLE			25	12
70	ADVA6763	GREEN	2-0	75	1/2 CIRCLE			20	12
71	ADVA6767	GREEN	2-0	90	1/2 CIRCLE			26	12

 Round Bodied

 Taper Cut

PL - Pledget

H - Heavy

Sr. No.	Double Arm Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
72	ADVA6817	WHITE	2-0	75	1/2 CIRCLE			17	12
73	ADVA6917	WHITE	2-0	90	1/2 CIRCLE			17	12
74	ADVA6917PL3	WHITE	2-0	90	1/2 CIRCLE			17	12
75	ADVA6917PL6	2-GREEN & 2-WHITE	2-0	75	1/2 CIRCLE			17	12
76	ADVA6917SPL3	WHITE	2-0	75	1/2 CIRCLE			17	12
77	ADVA6918PL6	2-GREEN & 2-WHITE	2-0	90	1/2 CIRCLE			25	12
78	ADVA6926	GREEN	2-0	90	1/2 CIRCLE			16	12
79	ADVA6927PL3	GREEN	2-0	75	1/2 CIRCLE			17	12
80	ADVA6937	GREEN	2-0	90	1/2 CIRCLE			17	12
81	ADVA6937S	GREEN	2-0	75	1/2 CIRCLE			17	12
82	ADVA6937PL3	GREEN	2-0	90	1/2 CIRCLE			17	12
83	ADVA6947PL3	WHITE	2-0	70	1/2 CIRCLE			17	12
84	ADVA6967	GREEN	2-0	90	1/2 CIRCLE			22	12
85	ADVA6977	GREEN	2-0	90	1/2 CIRCLE			26	12
86	ADVA6977PL3	GREEN	2-0	90	1/2 CIRCLE			26	12
87	ADVA6977PL6	GREEN	2-0	90	1/2 CIRCLE			26	12
88	ADVA6986PL6	WHITE	2-0	75	1/2 CIRCLE			25	12
89	ADVA6987	WHITE	2-0	90	1/2 CIRCLE			26	12
90	ADVA6987PL3	WHITE	2-0	90	1/2 CIRCLE			26	12
91	ADVA6987PL6	WHITE	2-0	90	1/2 CIRCLE			26	12
92	ADVA4B54PL3	4-GREEN & 4-WHITE	2-0	75	1/2 CIRCLE			16	12
93	ADVA10B52	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			17	6
94	ADVA10B55	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			17	6
95	ADVA10B55PL3	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			17	6
96	ADVA10B76	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			25	6
97	ADVA10B76PL6	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			25	6
98	ADVA10B77	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			25	6
99	ADVA10B77PL3	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			25	6
100	ADVA10B77PL6	5-GREEN & 5-WHITE	2-0	75	1/2 CIRCLE			26	6
101	ADVA655DA	GREEN	3-0	100	1/2 CIRCLE			25	12
102	ADVA6003	GREEN	3-0	90	1/2 CIRCLE			26	12
103	ADVA6006	GREEN	3-0	60	1/2 CIRCLE			16	12
104	ADVA6020	GREEN	3-0	90	3/8 CIRCLE			17	12
105	ADVA6026	GREEN	3-0	90	1/2 CIRCLE			18	12
106	ADVA6027	GREEN	3-0	75	1/2 CIRCLE			20	12
107	ADVA6169	GREEN	3-0	90	1/2 CIRCLE			16	12
108	ADVA6552	GREEN	3-0	100	1/2 CIRCLE	 LT	 LT	26	12
109	ADVA6759	GREEN	3-0	60	1/2 CIRCLE			17	12

 Round Bodied

 Taper Cut

PL - Pledget

H - Heavy

Sr. No.	Double Arm Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
110	ADVA6936	GREEN	3-0	90	1/2 CIRCLE			17	12
111	ADVA6936PL3	GREEN	3-0	90	1/2 CIRCLE			17	12
112	ADVA6976	GREEN	3-0	90	1/2 CIRCLE			26	12
113	ADVA18535	GREEN	3-0	45	3/8 CIRCLE			25	12
114	ADVA682	GREEN	4-0	76	1/2 CIRCLE			8	12
115	ADVA6021	WHITE	4-0	45	1/4 CIRCLE			8	12
116	ADVA6416	GREEN	4-0	90	1/2 CIRCLE			16	12
117	ADVA6835	GREEN	4-0	75	1/2 CIRCLE			17	12
118	ADVA6837	GREEN	4-0	75	1/2 CIRCLE			17	12
119	ADVA6891	GREEN	4-0	75	3/8 CIRCLE			13	12
120	ADVA6892PL3	GREEN	4-0	75	1/2 CIRCLE			13	12
121	ADVA6893PL3	WHITE	4-0	75	1/2 CIRCLE			13	12
122	ADVA6915	WHITE	4-0	60	3/8 CIRCLE			17	12
123	ADVA6935	GREEN	4-0	90	1/2 CIRCLE			17	12
124	ADVA18536	GREEN	4-0	45	3/8 CIRCLE			25	12
125	ADVA893	WHITE	5-0	45	1/4 CIRCLE			8	12
126	ADVA6005	GREEN	5-0	45	1/4 CIRCLE			8	12
127	ADVA6022	GREEN	5-0	45	1/4 CIRCLE			8	12
128	ADVA6545	GREEN	5-0	45	1/2 CIRCLE			8	12
129	ADVA6716	GREEN	5-0	75	1/2 CIRCLE			16	12
130	ADVA6757	GREEN	5-0	60	1/2 CIRCLE			17	12
131	ADVA6890	GREEN	5-0	75	3/8 CIRCLE			13	12
132	ADVA6890PL3	GREEN	5-0	75	3/8 CIRCLE			13	12
133	ADVA6597	GREEN	6-0	75	3/8 CIRCLE			11	12

⊙ Round Bodied

⊖ Taper Cut

△ Cutting

▽ Spatulated Micro Point

PL - Pledget

LT - Long Taper

Product Description


Composition	Polyethylene Terephthalate
Coating	Silicone
Sizes	Size 5 to 6-0
Colour	White (Undyed) & Green
Absorption Profile	Non-Absorbable
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Reverse Cutting, Taper Cut, Blunt Point, Oval Round Bodied, Diamond Point & Spatula Micro Point needles

Tensile Strength

ADVABOND Suture elicits a minimal initial inflammatory reaction in tissues, which is followed by a gradual encapsulation of the suture by fibrous connective tissue.¹

No significant change in the retention of tensile strength of the suture occurs during the lifetime of the implantation.¹

Knot pull tensile strength (KgF)²



REFERENCE _____

1 : Data on file with AMS - IFU ADVABOND

2 : Data on file with AMS - DV done on Polyester, size - 2-0, 17 mm Round bodied needle

Multifilament Material




Green



White (Undyed)

Premium Needles for consistent Superior Performance



Minimal tissue trauma enables faster recovery



P-21-22, 25-26 & 34-35, GIDC, Manjusar, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India
T : (91) 2667 671302 | amsltd.com | info@amsltd.com

Singapore Office : 111, Somerset Road, 07-08, 111 Somerset, Singapore - 238164.
T: (65) 636910 80 / 89 | E: info@amsltd.com | amsltd.com