



# ADVAMRYL

## Poliglecaprone - 25 (Undyed & Violet)

### Monofilament Synthetic Absorbable Surgical Suture

**Supreme handling performance for confidence  
in subcuticular skin closure**



**ADVATRAY™ - the confidence to  
suture with ease**



## Predictable strength for superior support

### Product Description

| Sr. No. | Code      | Suture |      |             | Needle     |                                                                                      |             |             | Units Per Box |
|---------|-----------|--------|------|-------------|------------|--------------------------------------------------------------------------------------|-------------|-------------|---------------|
|         |           | Colour | Size | Length (cm) | Curvature  | Type                                                                                 | Length (mm) | Speciality  |               |
| 1       | ADVA1047  | VIOLET | 2-0  | 150         | -          | -                                                                                    | -           | -           | 12            |
| 2       | ADVA1011  | UNDYED | 0    | 70          | 3/8 CIRCLE |    | 30          | -           | 12            |
| 3       | ADVA1014  | UNDYED | 0    | 90          | 1/2 CIRCLE |    | 48          | -           | 12            |
| 4       | ADVA1015  | VIOLET | 0    | 70          | 1/2 CIRCLE |    | 37          | -           | 12            |
| 5       | ADVA1026  | UNDYED | 0    | 70          | 1/2 CIRCLE |    | 40          | HEAVY       | 12            |
| 6       | ADVA1027  | VIOLET | 0    | 76          | 1/2 CIRCLE |    | 35          | -           | 12            |
| 7       | ADVA1029  | VIOLET | 0    | 70          | 1/2 CIRCLE |    | 37          | -           | 12            |
| 8       | ADVA1030  | VIOLET | 0    | 90          | 1/2 CIRCLE |    | 40          | HEAVY       | 12            |
| 9       | ADVA1041  | VIOLET | 0    | 70          | 1/2 CIRCLE |    | 40          | HEAVY       | 12            |
| 10      | ADVA1051  | VIOLET | 0    | 75          | 3/8 CIRCLE |    | 20          | -           | 12            |
| 11      | ADVA1052  | VIOLET | 0    | 70          | 1/2 CIRCLE |   | 30          | -           | 12            |
| 12      | ADVA1652  | UNDYED | 0    | 70          | 1/2 CIRCLE |  | 25          | -           | 12            |
| 13      | ADVA1642  | UNDYED | 0    | 70          | 1/2 CIRCLE |  | 30          | -           | 12            |
| 14      | ADVA1742  | VIOLET | 0    | 70          | 1/2 CIRCLE |  | 30          | -           | 12            |
| 15      | ADVA3442  | VIOLET | 0    | 70          | 1/2 CIRCLE |  | 31          | LONG TAPER  | 12            |
| 16      | ADVA3443  | VIOLET | 0    | 70          | 1/2 CIRCLE |  | 36          | LONG TAPER  | 12            |
| 17      | ADVA3601  | VIOLET | 0    | 70          | 1/2 CIRCLE |  | 36          | -           | 12            |
| 18      | ADVA1006  | UNDYED | 1    | 90          | 1/2 CIRCLE |  | 40          | -           | 12            |
| 19      | ADVA1016  | VIOLET | 1    | 70          | 1/2 CIRCLE |  | 37          | -           | 12            |
| 20      | ADVA1017  | UNDYED | 1    | 90          | 1/2 CIRCLE |  | 48          | -           | 12            |
| 21      | ADVA1031  | VIOLET | 1    | 70          | 1/2 CIRCLE |  | 40          | HEAVY       | 12            |
| 22      | ADVA1042  | VIOLET | 1    | 90          | 1/2 CIRCLE |  | 40          | HEAVY       | 12            |
| 23      | ADVA1043  | VIOLET | 1    | 90          | 1/2 CIRCLE |  | 50          | BLUNT POINT | 12            |
| 24      | ADVA3457  | VIOLET | 1    | 70          | 1/2 CIRCLE |  | 36          | LONG TAPER  | 12            |
| 25      | ADVA3755  | VIOLET | 1    | 70          | 1/2 CIRCLE |  | 31          | LONG TAPER  | 12            |
| 26      | ADVA3759U | UNDYED | 1    | 70          | 1/2 CIRCLE |  | 40          | -           | 12            |
| 27      | ADVA945P  | UNDYED | 2-0  | 90          | 1/2 CIRCLE |  | 36          | -           | 12            |
| 28      | ADVA1005  | UNDYED | 2-0  | 70          | 3/8 CIRCLE |  | 30          | -           | 12            |
| 29      | ADVA1018  | VIOLET | 2-0  | 70          | 3/8 CIRCLE |  | 24          | -           | 12            |
| 30      | ADVA1019  | VIOLET | 2-0  | 70          | 1/2 CIRCLE |  | 26          | LONG TAPER  | 12            |
| 31      | ADVA1020  | VIOLET | 2-0  | 70          | 1/2 CIRCLE |  | 37          | -           | 12            |
| 32      | ADVA1044  | VIOLET | 2-0  | 70          | 1/2 CIRCLE |  | 35          | -           | 12            |
| 33      | ADVA1053  | VIOLET | 2-0  | 70          | 1/2 CIRCLE |  | 30          | -           | 12            |

 Round Bodied

 Oval Round Bodied

 Reverse Cutting

 Taper Cut

| Sr. No. | Code       | Suture |      |             | Needle     |                                                                                      |             |            | Units Per Box |
|---------|------------|--------|------|-------------|------------|--------------------------------------------------------------------------------------|-------------|------------|---------------|
|         |            | Colour | Size | Length (cm) | Curvature  | Type                                                                                 | Length (mm) | Speciality |               |
| 34      | ADVA1075   | VIOLET | 2-0  | 75          | 1/2 CIRCLE |    | 26          | LONG TAPER | 12            |
| 35      | ADVA1603   | UNDYED | 2-0  | 70          | 1/2 CIRCLE |    | 25          | -          | 12            |
| 36      | ADVA1602   | UNDYED | 2-0  | 70          | 5/8 CIRCLE |    | 26          | -          | 12            |
| 37      | ADVA1641   | UNDYED | 2-0  | 70          | 1/2 CIRCLE |    | 30          | -          | 12            |
| 38      | ADVA1665   | UNDYED | 2-0  | 70          | 1/2 CIRCLE |    | 26          | -          | 12            |
| 39      | ADVA1666   | UNDYED | 2-0  | 70          | 1/2 CIRCLE |    | 36          | -          | 12            |
| 40      | ADVA1666S  | UNDYED | 2-0  | 45          | 1/2 CIRCLE |    | 36          | -          | 12            |
| 41      | ADVA1741   | VIOLET | 2-0  | 70          | 1/2 CIRCLE |    | 30          | -          | 12            |
| 42      | ADVA3200   | VIOLET | 2-0  | 70          | 1/2 CIRCLE |    | 31          | LONG TAPER | 12            |
| 43      | ADVA3327   | UNDYD  | 2-0  | 70          | 3/8 CIRCLE |    | 26          | -          | 12            |
| 44      | ADVA3431   | UNDYD  | 2-0  | 70          | 3/8 CIRCLE |    | 26          | -          | 12            |
| 45      | ADVA3441   | VIOLET | 2-0  | 70          | 1/2 CIRCLE |    | 31          | -          | 12            |
| 46      | ADVA3651   | UNDYED | 2-0  | 70          | STRAIGHT   |    | 60          | -          | 12            |
| 47      | ADVA3663   | VIOLET | 2-0  | 70          | 1/2 CIRCLE |    | 19          | -          | 12            |
| 48      | ADVA1001   | UNDYED | 3-0  | 70          | 1/2 CIRCLE |    | 20          | -          | 12            |
| 49      | ADVA1004   | UNDYED | 3-0  | 70          | 3/8 CIRCLE |    | 24          | -          | 12            |
| 50      | ADVA1008   | UNDYED | 3-0  | 90          | 3/8 CIRCLE |    | 19          | -          | 12            |
| 51      | ADVA1010   | UNDYED | 3-0  | 75          | 3/8 CIRCLE |   | 26          | -          | 12            |
| 52      | ADVA1021SA | UNDYED | 3-0  | 70          | 1/2 CIRCLE |  | 26          | -          | 12            |
| 53      | ADVA1022   | VIOLET | 3-0  | 70          | 3/8 CIRCLE |  | 24          | -          | 12            |
| 54      | ADVA1023   | VIOLET | 3-0  | 70          | 1/2 CIRCLE |  | 30          | -          | 12            |
| 55      | ADVA1028   | UNDYED | 3-0  | 70          | 1/2 CIRCLE |  | 23          | -          | 12            |
| 56      | ADVA1032   | UNDYED | 3-0  | 70          | 3/8 CIRCLE |  | 24          | -          | 12            |
| 57      | ADVA1033   | VIOLET | 3-0  | 70          | 1/2 CIRCLE |  | 26          | LONG TAPER | 12            |
| 58      | ADVA1035   | VIOLET | 3-0  | 70          | 1/2 CIRCLE |  | 22          | LONG TAPER | 12            |
| 59      | ADVA1046   | UNDYED | 3-0  | 70          | 1/2 CIRCLE |  | 26          | -          | 12            |
| 60      | ADVA1049   | VIOLET | 3-0  | 90          | 3/8 CIRCLE |  | 19          | PRIME      | 12            |
| 61      | ADVA1050   | VIOLET | 3-0  | 75          | STRAIGHT   |  | 60          | -          | 12            |
| 62      | ADVA1054   | VIOLET | 3-0  | 70          | 3/8 CIRCLE |  | 20          | PRIME      | 12            |
| 63      | ADVA1056   | UNDYED | 3-0  | 45          | 3/8 CIRCLE |  | 19          | PRIME      | 12            |
| 64      | ADVA1064   | UNDYED | 3-0  | 45          | 3/8 CIRCLE |  | 19          | -          | 12            |
| 65      | ADVA1065   | VIOLET | 3-0  | 70          | 3/8 CIRCLE |  | 19          | -          | 12            |
| 66      | ADVA1067   | UNDYED | 3-0  | 70          | 3/8 CIRCLE |  | 16          | -          | 12            |
| 67      | ADVA1070   | VIOLET | 3-0  | 75          | 3/8 CIRCLE |  | 19          | -          | 12            |
| 68      | ADVA1077   | VIOLET | 3-0  | 75          | 3/8 CIRCLE |  | 26          | -          | 12            |
| 69      | ADVA1079   | UNDYED | 3-0  | 45          | 3/8 CIRCLE |  | 16          | -          | 12            |
| 70      | ADVA1326   | UNDYED | 3-0  | 70          | 3/8 CIRCLE |  | 25          | PRIME      | 12            |
| 71      | ADVA1326S  | UNDYED | 3-0  | 45          | 3/8 CIRCLE |  | 25          | PRIME      | 12            |
| 72      | ADVA1326MR | VIOLET | 3-0  | 70          | 3/8 CIRCLE |  | 25          | PRIME      | 12            |

 Round Bodied

 Oval Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut



| Sr. No. | Code       | Suture |      |             | Needle     |                                                                                      |             |            | Units Per Box |
|---------|------------|--------|------|-------------|------------|--------------------------------------------------------------------------------------|-------------|------------|---------------|
|         |            | Colour | Size | Length (cm) | Curvature  | Type                                                                                 | Length (mm) | Speciality |               |
| 73      | ADVA1663   | UNDYED | 3-0  | 70          | 1/2 CIRCLE |    | 22          | -          | 12            |
| 74      | ADVA1664   | UNDYED | 3-0  | 70          | 1/2 CIRCLE |    | 26          | -          | 12            |
| 75      | ADVA1737   | VIOLET | 3-0  | 70          | 1/2 CIRCLE |    | 20          | -          | 12            |
| 76      | ADVA3202   | UNDYED | 3-0  | 70          | 3/8 CIRCLE |    | 19          | -          | 12            |
| 77      | ADVA3207   | UNDYED | 3-0  | 45          | 3/8 CIRCLE |    | 19          | PRIME      | 12            |
| 78      | ADVA3326M  | UNDYED | 3-0  | 70          | 3/8 CIRCLE |    | 26          | -          | 12            |
| 79      | ADVA3326MS | VIOLET | 3-0  | 45          | 3/8 CIRCLE |    | 25          | -          | 12            |
| 80      | ADVA3650   | UNDYED | 3-0  | 70          | STRAIGHT   |    | 60          | -          | 12            |
| 81      | ADVA3665   | VIOLET | 3-0  | 75          | 1/2 CIRCLE |    | 31          | LONG TAPER | 12            |
| 82      | ADVA1003   | UNDYED | 4-0  | 70          | 3/8 CIRCLE |    | 24          | -          | 12            |
| 83      | ADVA1007   | UNDYED | 4-0  | 75          | 1/2 CIRCLE |    | 17          | -          | 12            |
| 84      | ADVA1009   | UNDYED | 4-0  | 45          | 3/8 CIRCLE |    | 19          | PRIME      | 12            |
| 85      | ADVA1024   | VIOLET | 4-0  | 70          | 3/8 CIRCLE |    | 19          | -          | 12            |
| 86      | ADVA1025   | VIOLET | 4-0  | 70          | 1/2 CIRCLE |    | 17          | -          | 12            |
| 87      | ADVA1034   | VIOLET | 4-0  | 70          | 1/2 CIRCLE |    | 22          | LONG TAPER | 12            |
| 88      | ADVA1036   | UNDYED | 4-0  | 70          | 3/8 CIRCLE |    | 19          | -          | 12            |
| 89      | ADVA1038   | VIOLET | 4-0  | 75          | 1/2 CIRCLE |    | 20          | LONG TAPER | 12            |
| 90      | ADVA1048   | UNDYED | 4-0  | 70          | 3/8 CIRCLE |   | 19          | PRIME      | 12            |
| 91      | ADVA1058   | VIOLET | 4-0  | 45          | 3/8 CIRCLE |  | 19          | -          | 12            |
| 92      | ADVA1066   | UNDYED | 4-0  | 70          | 3/8 CIRCLE |  | 12          | -          | 12            |
| 93      | ADVA1069   | UNDYED | 4-0  | 45          | 3/8 CIRCLE |  | 13          | -          | 12            |
| 94      | ADVA1073   | VIOLET | 4-0  | 50          | 3/8 CIRCLE |  | 19          | ADVAPASS   | 12            |
| 95      | ADVA1076   | UNDYED | 4-0  | 45          | 3/8 CIRCLE |  | 12          | -          | 12            |
| 96      | ADVA1078   | UNDYED | 4-0  | 75          | 3/8 CIRCLE |  | 20          | -          | 12            |
| 97      | ADVA1134   | UNDYED | 4-0  | 70          | 1/2 CIRCLE |  | 20          | -          | 12            |
| 98      | ADVA1205   | UNDYED | 4-0  | 70          | 3/8 CIRCLE |  | 16          | PRIME      | 12            |
| 99      | ADVA1205MR | VIOLET | 4-0  | 70          | 3/8 CIRCLE |  | 16          | PRIME      | 12            |
| 100     | ADVA1648   | UNDYED | 4-0  | 70          | 3/8 CIRCLE |  | 16          | LONG TAPER | 12            |
| 101     | ADVA1748   | VIOLET | 4-0  | 70          | 3/8 CIRCLE |  | 16          | LONG TAPER | 12            |
| 102     | ADVA3201   | UNDYED | 4-0  | 45          | 3/8 CIRCLE |  | 19          | -          | 12            |
| 103     | ADVA3205   | UNDYED | 4-0  | 45          | 3/8 CIRCLE |  | 16          | PRIME      | 12            |
| 104     | ADVA845G   | UNDYED | 4-0  | 45          | 3/8 CIRCLE |  | 16          | PRIME      | 12            |
| 105     | ADVA1002   | UNDYED | 5-0  | 45          | 3/8 CIRCLE |  | 19          | PRIME      | 12            |
| 106     | ADVA1037   | VIOLET | 5-0  | 70          | 1/2 CIRCLE |  | 17          | -          | 12            |
| 107     | ADVA1039   | VIOLET | 5-0  | 75          | 3/8 CIRCLE |  | 13          | LONG TAPER | 12            |
| 108     | ADVA1055   | VIOLET | 5-0  | 40          | CURVED     |  | 20          | -          | 12            |
| 109     | ADVA1057   | UNDYED | 5-0  | 70          | 3/8 CIRCLE |  | 12          | PRIME      | 12            |
| 110     | ADVA1059   | VIOLET | 5-0  | 45          | 3/8 CIRCLE |  | 16          | -          | 12            |
| 111     | ADVA1060   | VIOLET | 5-0  | 45          | 3/8 CIRCLE |  | 16          | -          | 12            |

|                                                                                                  |                                                                                                       |                                                                                                     |                                                                                               |                                                                                                 |
|--------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
|  Round Bodied |  Oval Round Bodied |  Reverse Cutting |  Cutting |  Taper Cut |
|--------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|

| Sr. No. | Code     | Suture |      |             | Needle     |                                                                                     |             |            | Units Per Box |
|---------|----------|--------|------|-------------|------------|-------------------------------------------------------------------------------------|-------------|------------|---------------|
|         |          | Colour | Size | Length (cm) | Curvature  | Type                                                                                | Length (mm) | Speciality |               |
| 112     | ADVA1062 | UNDYED | 5-0  | 70          | 3/8 CIRCLE |   | 16          | -          | 12            |
| 113     | ADVA1068 | VIOLET | 5-0  | 70          | 3/8 CIRCLE |   | 18          | -          | 12            |
| 114     | ADVA1072 | VIOLET | 5-0  | 45          | 3/8 CIRCLE |   | 12          | -          | 12            |
| 115     | ADVA1074 | VIOLET | 5-0  | 50          | 3/8 CIRCLE |   | 16          | -          | 12            |
| 116     | ADVA3203 | UNDYED | 5-0  | 45          | 3/8 CIRCLE |   | 13          | PRIME      | 12            |
| 117     | ADVA3204 | UNDYED | 5-0  | 45          | 3/8 CIRCLE |   | 16          | PRIME      | 12            |
| 118     | ADVA3209 | UNDYED | 5-0  | 70          | 3/8 CIRCLE |   | 19          | -          | 12            |
| 119     | ADVA493G | UNDYED | 5-0  | 45          | 3/8 CIRCLE |   | 13          | PRIME      | 12            |
| 120     | ADVA844G | UNDYED | 5-0  | 45          | 3/8 CIRCLE |   | 16          | PRIME      | 12            |
| 121     | ADVA1040 | VIOLET | 6-0  | 75          | 1/2 CIRCLE |   | 13          | LONG TAPER | 12            |
| 122     | ADVA1061 | UNDYED | 6-0  | 70          | 3/8 CIRCLE |   | 12          | PRIME      | 12            |
| 123     | ADVA1063 | UNDYED | 6-0  | 45          | 3/8 CIRCLE |   | 13          | -          | 12            |
| 124     | ADVA1071 | UNDYED | 6-0  | 50          | 3/8 CIRCLE |   | 12          | PRIME      | 12            |
| 125     | ADVA3214 | UNDYED | 6-0  | 45          | 3/8 CIRCLE |   | 11          | PRIME      | 12            |
| 126     | ADVA3224 | VIOLET | 6-0  | 45          | 1/2 CIRCLE |  | 13          | -          | 12            |

| Sr. No. | Code     | Suture |      |             | Needle                  |                                                                                      |                                                                                       |             | Units Per Box |
|---------|----------|--------|------|-------------|-------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|-------------|---------------|
|         |          | Colour | Size | Length (cm) | Curvature               | Type 1                                                                               | Type 2                                                                                | Length (mm) |               |
| 127     | ADVA1012 | VIOLET | 0    | 75          | 1/2 CIRCLE              |  |  | 30          | 12            |
| 128     | ADVA1021 | VIOLET | 3-0  | 70          | 1/2 CIRCLE & 3/8 CIRCLE |  |  | 26          | 12            |

 Round Bodied

 Reverse Cutting

 Cutting

### Product Description


|                     |                                                                                                                                                               |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Composition         | Glycolide and ̵-Caprolactone                                                                                                                                  |
| Coating             | None                                                                                                                                                          |
| Size                | Size 1 to 6-0                                                                                                                                                 |
| Colour              | Undyed and Violet                                                                                                                                             |
| Absorption Profile  | Essentially complete up to 90 days                                                                                                                            |
| Needle Availability | A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Trocar, Reverse Cutting, Taper Cut, Blunt Point, Oval Round Bodied & Diamond Point needles |

### Tensile Strength

| Days of Implantation     | Tensile Strength Retention |        |
|--------------------------|----------------------------|--------|
|                          | DYED                       | UNDYED |
| 7 Days (1 to 5-0)        | 60%                        | 55%    |
| 7 Days (6-0)             | 60%                        | 50%    |
| 14 Days (1 to 3-0 & 5-0) | 40%                        | 20%    |
| 14 Days (4-0 & 6-0)      | 40%                        | 15%    |

66% more knot pull tensile strength Vs. recommended limit ensures excellent support throughout the critical wound healing period


### Knot pull tensile strength (KgF)<sup>1</sup>




REFERENCE \_\_\_\_\_

1 : Data on file with AMS - Bench test done on Poliglecaprone - 25, size 3-0,

### Monofilament Material



Violet



Undyed


WEEKS



Provides excellent support for a range of critical wound healing period to secure closure in a variety of tissues, procedures and patients.

# ADVA TRAY™ - the confidence to suture with ease

**ADVA TRAY™** a unique dispensing system designed to deliver ease of suture dispensing with minimal memory in the suture strand



## Premium Needles for consistent Superior Performance

**PREMIUM NEEDLES**  
ensure  
**UNSURPASSED SUTURING EXPERIENCE**

|                                                                |                                                                                                                                                         |
|----------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>DURABLE COAT</b><br>RETAINS SHARPNESS...<br>PASS AFTER PASS | Maintains needle sharpness over the duration of procedure, resulting into consistent penetration pass after pass thus less potential for tissue trauma. |
| <b>RIBBED BODY</b><br>FOR SUPERIOR CONTROL                     | For superior grasping and stability of needle in needle forceps.                                                                                        |
| <b>CHAMFERING</b><br>MINIMUM SUTURE DETACHMENT                 | Diminishes chances for suture fracturing at the swage by rounding off the sharp edges of the barrel in the swage area of the needle.                    |
| <b>ELECTROPOLISHING</b><br>SMOOTH & BRIGHT                     | Makes needle surface free from micro asperities, thus creates ultra sharp edges for lesser tissue drag.                                                 |

Minimal tissue trauma enables faster recovery



P-21-22, 25-26 & 34-35, GIDC, Manjusar, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India  
T : (91) 2667 671302 | [amsltd.com](http://amsltd.com) | [info@amsltd.com](mailto:info@amsltd.com)

Singapore Office : 111, Somerset Road, 07-08, 111 Somerset, Singapore - 238164.  
T: (65) 636910 80 / 89 | E: [info@amsltd.com](mailto:info@amsltd.com) | [amsltd.com](http://amsltd.com)