



ADVALON Polyamide

Monofilament Synthetic Non-Absorbable Surgical Suture

**Ideal selection for
confidence in skin closure**



Appropriate choice for dermal closure and retention suturing

Product Description

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVA905	BLACK	1	2X76	-	-	-	-	12
2	ADVA739	BLACK	0	150	1/2 CIRCLE	 LP	31	-	12
3	ADVA3004	BLACK	0	90	3/8 CIRCLE		40	-	12
4	ADVA3011	BLACK	0	90	1/2 CIRCLE		36	-	12
5	ADVA3035	BLACK	0	75	3/8 CIRCLE		24	-	12
6	ADVA3036	BLACK	0	75	3/8 CIRCLE		30	-	12
7	ADVA3053	BLACK	0	150	1/2 CIRCLE	 LP	50	HEAVY	12
8	ADVA3059	BLACK	0	75	3/8 CIRCLE		30	-	12
9	ADVA3063	BLACK	0	150	1/2 CIRCLE		40	-	12
10	ADVA3064	BLACK	0	150	3/8 CIRCLE		40	-	12
11	ADVA3118	BLACK	0	100	1/2 CIRCLE		40	-	12
12	ADVA3119	BLACK	0	90	3/8 CIRCLE		24	-	12
13	ADVA3120	BLACK	0	90	3/8 CIRCLE		30	-	12
14	ADVA3121	BLACK	0	45	3/8 CIRCLE		30	-	12
15	ADVA3122	BLACK	0	100	3/8 CIRCLE		60	-	12
16	ADVA3123	BLACK	0	75	STRAIGHT		75	-	12
17	ADVA3227	BLACK	0	70	3/8 CIRCLE		36	-	12
18	ADVA3256	BLACK	0	75	3/8 CIRCLE		39	-	12
19	ADVA3304	BLACK	0	75	1/2 CIRCLE		38	-	12
20	ADVA3334	BLACK	0	100	1/2 CIRCLE		40	HEAVY	12
21	ADVA3337	BLACK	0	70	3/8 CIRCLE		45	-	12
22	ADVA3340	BLACK	0	150	1/2 CIRCLE	 LP	40	-	12
23	ADVA3346	BLACK	0	100	1/2 CIRCLE		40	HEAVY	12
24	ADVA3389	BLACK	0	90	STRAIGHT		60	-	12
25	ADVA3741	BLACK	0	150	1/2 CIRCLE	 LP	48	HEAVY	12
26	ADVA738	BLACK	1	100	1/2 CIRCLE		48	-	12
27	ADVA743	BLACK	1	150	1/2 CIRCLE	 LP	36	-	12
28	ADVA3006	BLACK	1	75	3/8 CIRCLE		90	-	12
29	ADVA3012	BLACK	1	90	1/2 CIRCLE		48	-	12
30	ADVA3032	BLACK	1	100	1/2 CIRCLE		90	-	12
31	ADVA3033	BLACK	1	150	1/2 CIRCLE	 LP	50	HEAVY	12
32	ADVA3037	BLACK	1	75	3/8 CIRCLE		30	-	12
33	ADVA3038	BLACK	1	75	1/2 CIRCLE		30	-	12
34	ADVA3065	BLACK	1	90	1/2 CIRCLE		40	-	12
35	ADVA3075	BLACK	1	100	1/2 CIRCLE	 LP	50	HEAVY	12

 Round Bodied

 Reverse Cutting

 Cutting

 LP - Loop

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
36	ADVA3082	BLACK	1	90	1/2 CIRCLE		35	-	12
37	ADVA3086	BLACK	1	90	1/2 CIRCLE		40	-	12
38	ADVA3145	BLACK	1	90	1/2 CIRCLE		45	-	12
39	ADVA3114	BLACK	1	100	3/8 CIRCLE		90	-	12
40	ADVA3117	BLACK	1	75	1/2 CIRCLE		37	-	12
41	ADVA3124	BLACK	1	75	3/8 CIRCLE		39	-	12
42	ADVA3125	BLACK	1	100	1/2 CIRCLE		30	-	12
43	ADVA3126	BLACK	1	100	1/2 CIRCLE		43	HEAVY	12
44	ADVA3127	BLACK	1	75	1/2 CIRCLE		48	-	12
45	ADVA3128	BLACK	1	45	3/8 CIRCLE		30	-	12
46	ADVA3129	BLACK	1	90	3/8 CIRCLE		30	-	12
47	ADVA3130	BLACK	1	90	STRAIGHT		51	-	12
48	ADVA3170	BLACK	1	90	1/2 CIRCLE		40	-	12
49	ADVA3171	BLACK	1	75	3/8 CIRCLE		30	-	12
50	ADVA3175	BLACK	1	150	3/8 CIRCLE		75	-	12
51	ADVA3305	BLACK	1	75	1/2 CIRCLE		38	-	12
52	ADVA3335	BLACK	1	100	1/2 CIRCLE		40	HEAVY	12
53	ADVA3338	BLACK	1	100	1/2 CIRCLE		45	HEAVY	12
54	ADVA3347	BLACK	1	100	1/2 CIRCLE		40	HEAVY	12
55	ADVA3348	BLACK	1	150	1/2 CIRCLE	 LP	48	HEAVY	12
56	ADVA3348SP	BLACK	1	150	1/2 CIRCLE		40	HEAVY	12
57	ADVA3362	BLACK	1	100	1/2 CIRCLE		50	HEAVY	12
58	ADVA3397	BLACK	1	70	3/8 CIRCLE		90	-	12
59	ADVA3747	BLACK	1	100	1/2 CIRCLE		40	-	12
60	ADVA3057	BLACK	2	75	1/2 CIRCLE		40	-	12
61	ADVA3076	BLACK	2	90	1/2 CIRCLE		37	-	12
62	ADVA3087	BLACK	2	90	1/2 CIRCLE		40	-	12
63	ADVA3131	BLACK	2	100	3/8 CIRCLE		60	-	12
64	ADVA3306	BLACK	2	75	1/2 CIRCLE		38	-	12
65	ADVA568	BLACK	2-0	100	1/2 CIRCLE		31	-	12
66	ADVA0786	BLUE	2-0	45	3/8 CIRCLE		40	-	12
67	ADVA3002	BLACK	2-0	90	3/8 CIRCLE		40	-	12
68	ADVA3002L	BLACK	2-0	100	3/8 CIRCLE		40	-	12
69	ADVA3007	BLACK	2-0	100	3/8 CIRCLE		45	-	12
70	ADVA3010	BLACK	2-0	90	3/8 CIRCLE		24	-	12
71	ADVA3013	BLACK	2-0	90	3/8 CIRCLE		36	-	12
72	ADVA3014	BLACK	2-0	90	3/8 CIRCLE		19	-	12
73	ADVA3017	BLACK	2-0	100	STRAIGHT		55	-	12
74	ADVA3021	BLACK	2-0	100	1/2 CIRCLE		30	-	12
75	ADVA3027	BLACK	2-0	100	1/2 CIRCLE		45	-	12
76	ADVA3027SP	BLACK	2-0	76	1/2 CIRCLE		50	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
77	ADVA3039	BLACK	2-0	75	3/8 CIRCLE		24	-	12
78	ADVA3044	BLACK	2-0	75	3/8 CIRCLE		24	-	12
79	ADVA3054	BLACK	2-0	45	1/2 CIRCLE		40	-	12
80	ADVA3056	BLACK	2-0	75	3/8 CIRCLE		31	-	12
81	ADVA3069	BLACK	2-0	90	3/8 CIRCLE		60	PRIME	12
82	ADVA3073	BLACK	2-0	75	1/2 CIRCLE		37	-	12
83	ADVA3074	BLACK	2-0	75	3/8 CIRCLE		39	-	12
84	ADVA3077	BLACK	2-0	76	3/8 CIRCLE		26	-	12
85	ADVA3088	BLACK	2-0	75	STRAIGHT		60	-	12
86	ADVA3096	BLACK	2-0	45	3/8 CIRCLE		24	-	12
87	ADVA3097	BLACK	2-0	45	3/8 CIRCLE		19	PRIME	12
88	ADVA3105	BLACK	2-0	45	3/8 CIRCLE		26	-	12
89	ADVA3107	BLACK	2-0	45	3/8 CIRCLE		26	-	12
90	ADVA3112	BLACK	2-0	45	3/8 CIRCLE		40	-	12
91	ADVA3132	BLACK	2-0	45	3/8 CIRCLE		30	-	12
92	ADVA3133	BLACK	2-0	90	STRAIGHT		51	-	12
93	ADVA3134	BLACK	2-0	75	STRAIGHT		60	-	12
94	ADVA3135	BLACK	2-0	45	3/8 CIRCLE		34	MICRO POINT	12
95	ADVA3136	BLACK	2-0	75	3/8 CIRCLE		24	MICRO POINT	12
96	ADVA3168	BLACK	2-0	90	3/8 CIRCLE		30	-	12
97	ADVA3169	BLACK	2-0	100	3/8 CIRCLE		60	-	12
98	ADVA3172	BLACK	2-0	90	3/8 CIRCLE		26	-	12
99	ADVA3173	BLACK	2-0	75	1/2 CIRCLE		26	-	12
100	ADVA3174	BLACK	2-0	90	1/2 CIRCLE		26	-	12
101	ADVA3176	BLACK	2-0	90	1/2 CIRCLE		30	-	12
102	ADVA3177	BLACK	2-0	90	1/2 CIRCLE		37	-	12
103	ADVA3230	BLACK	2-0	75	3/8 CIRCLE		30	-	12
104	ADVA3258	BLACK	2-0	100	3/8 CIRCLE		58	-	12
105	ADVA3275	BLACK	2-0	75	3/8 CIRCLE		26	-	12
106	ADVA3333	BLACK	2-0	70	3/8 CIRCLE		40	-	12
107	ADVA3333L	BLACK	2-0	90	3/8 CIRCLE		40	-	12
108	ADVA3336	BLACK	2-0	70	3/8 CIRCLE		45	-	12
109	ADVA3336DK	BLACK	2-0	110	3/8 CIRCLE		45	-	12
110	ADVA3336L	BLACK	2-0	90	3/8 CIRCLE		45	-	12
111	ADVA3336S	BLACK	2-0	35	3/8 CIRCLE		45	-	12
112	ADVA3336SP	BLACK	2-0	75	1/2 CIRCLE		30	-	12
113	ADVA3390	BLACK	2-0	90	STRAIGHT		60	-	12
114	ADVA3786	BLACK	2-0	45	3/8 CIRCLE		40	-	12
115	ADVA320	BLACK	3-0	45	3/8 CIRCLE		26	-	12
116	ADVA0663	BLACK	3-0	45	3/8 CIRCLE		24	-	12
117	ADVA3005	BLACK	3-0	70	3/8 CIRCLE		26	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
118	ADVA3009	BLACK	3-0	70	3/8 CIRCLE		19	-	12
119	ADVA3015	BLACK	3-0	90	3/8 CIRCLE		36	-	12
120	ADVA3028	BLACK	3-0	45	1/2 CIRCLE		40	-	12
121	ADVA3029	BLACK	3-0	45	1/2 CIRCLE		26	-	12
122	ADVA3034	BLACK	3-0	75	3/8 CIRCLE		24	-	12
123	ADVA3041	BLACK	3-0	75	3/8 CIRCLE		30	-	12
124	ADVA3042	BLACK	3-0	75	1/2 CIRCLE		30	-	12
125	ADVA3050	BLACK	3-0	75	3/8 CIRCLE		25	-	12
126	ADVA3068	BLACK	3-0	75	3/8 CIRCLE		26	PRIME	12
127	ADVA3071	BLACK	3-0	90	1/2 CIRCLE		24	-	12
128	ADVA3072	BLACK	3-0	40	3/8 CIRCLE		35	-	12
129	ADVA3078	BLACK	3-0	76	3/8 CIRCLE		19	-	12
130	ADVA3079	BLACK	3-0	76	3/8 CIRCLE		24	-	12
131	ADVA3083	BLACK	3-0	90	3/8 CIRCLE		24	-	12
132	ADVA3085	BLACK	3-0	90	1/2 CIRCLE		40	-	12
133	ADVA3089	BLACK	3-0	45	3/8 CIRCLE		20	-	12
134	ADVA3090	BLACK	3-0	75	STRAIGHT		60	-	12
135	ADVA3094	BLACK	3-0	45	3/8 CIRCLE		24	-	12
136	ADVA3101	BLACK	3-0	45	3/8 CIRCLE		16	PRIME	12
137	ADVA3104	BLACK	3-0	45	3/8 CIRCLE		18	-	12
138	ADVA3106	BLACK	3-0	45	3/8 CIRCLE		24	-	12
139	ADVA3137	BLACK	3-0	45	3/8 CIRCLE		19	-	12
140	ADVA3138	BLACK	3-0	75	3/8 CIRCLE		19	-	12
141	ADVA3139	BLACK	3-0	75	3/8 CIRCLE		39	-	12
142	ADVA3140	BLACK	3-0	90	3/8 CIRCLE		19	-	12
143	ADVA3141	BLACK	3-0	75	3/8 CIRCLE		16	-	12
144	ADVA3142	BLACK	3-0	90	3/8 CIRCLE		19	-	12
145	ADVA3143	BLACK	3-0	45	3/8 CIRCLE		30	-	12
146	ADVA3144	BLACK	3-0	90	3/8 CIRCLE		30	-	12
147	ADVA3146	BLACK	3-0	75	3/8 CIRCLE		39	-	12
148	ADVA3147	BLACK	3-0	90	STRAIGHT		51	-	12
149	ADVA3148	BLACK	3-0	90	STRAIGHT		60	-	12
150	ADVA3149	BLACK	3-0	45	3/8 CIRCLE		19	MICRO POINT	12
151	ADVA3150	BLACK	3-0	75	3/8 CIRCLE		19	MICRO POINT	12
152	ADVA3151	BLACK	3-0	45	3/8 CIRCLE		24	MICRO POINT	12
153	ADVA3152	BLACK	3-0	75	3/8 CIRCLE		24	MICRO POINT	12
154	ADVA3153	BLACK	3-0	45	3/8 CIRCLE		19	MICRO POINT	12
155	ADVA3178	BLACK	3-0	90	1/2 CIRCLE		21	-	12
156	ADVA3179	BLACK	3-0	90	1/2 CIRCLE		26	-	12
157	ADVA3180	BLACK	3-0	90	1/2 CIRCLE		30	-	12
158	ADVA3221	BLUE	3-0	90	3/8 CIRCLE		24	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
159	ADVA3243	BLACK	3-0	75	3/8 CIRCLE		24	-	12
160	ADVA3321	BLACK	3-0	70	3/8 CIRCLE		16	PRIME	12
161	ADVA3321SP	BLACK	3-0	70	3/8 CIRCLE		30	-	12
162	ADVA3328	BLACK	3-0	70	3/8 CIRCLE		26	-	12
163	ADVA3328G	BLACK	3-0	70	3/8 CIRCLE		26	-	12
164	ADVA3328L	BLACK	3-0	90	3/8 CIRCLE		26	-	12
165	ADVA3328S	BLACK	3-0	35	3/8 CIRCLE		26	-	12
166	ADVA3328SG	BLACK	3-0	35	3/8 CIRCLE		26	-	12
167	ADVA3332	BLACK	3-0	70	3/8 CIRCLE		16	-	12
168	ADVA9136	BLACK	3-0	90	1/2 CIRCLE		26	-	12
169	ADVAN662	BLACK	4-0	75	3/8 CIRCLE		19	-	12
170	ADVAN665	BLACK	4-0	45	3/8 CIRCLE		19	-	12
171	ADVAN667	BLACK	4-0	45	3/8 CIRCLE		12	-	12
172	ADVA1162	BLACK	4-0	45	1/2 CIRCLE		16	-	12
173	ADVA3001	BLACK	4-0	90	1/2 CIRCLE		16	-	12
174	ADVA3030	BLACK	4-0	45	1/2 CIRCLE		19	-	12
175	ADVA3043	BLACK	4-0	75	3/8 CIRCLE		24	-	12
176	ADVA3051	BLACK	4-0	75	3/8 CIRCLE		19	-	12
177	ADVA3055	BLACK	4-0	75	1/2 CIRCLE		16	-	12
178	ADVA3061	BLACK	4-0	90	3/8 CIRCLE		16	-	12
179	ADVA3062	BLACK	4-0	90	3/8 CIRCLE		20	-	12
180	ADVA3080	BLACK	4-0	76	3/8 CIRCLE		19	-	12
181	ADVA3084	BLACK	4-0	70	3/8 CIRCLE		16	-	12
182	ADVA3091	BLACK	4-0	45	3/8 CIRCLE		15	-	12
183	ADVA3095	BLACK	4-0	45	3/8 CIRCLE		16	-	12
184	ADVA3098	BLACK	4-0	75	3/8 CIRCLE		16	-	12
185	ADVA3099	BLACK	4-0	45	3/8 CIRCLE		24	-	12
186	ADVA3100	BLACK	4-0	45	3/8 CIRCLE		12	-	12
187	ADVA3108	BLACK	4-0	75	3/8 CIRCLE		30	-	12
188	ADVA3110	BLACK	4-0	45	3/8 CIRCLE		19	-	12
189	ADVA3116	BLACK	4-0	75	3/8 CIRCLE		20	SLIM BLADE	12
190	ADVA3154	BLACK	4-0	90	1/2 CIRCLE		21	-	12
191	ADVA3155	BLACK	4-0	90	3/8 CIRCLE		19	-	12
192	ADVA3156	BLACK	4-0	90	STRAIGHT		51	-	12
193	ADVA3157	BLACK	4-0	45	3/8 CIRCLE		13	MICRO POINT	12
194	ADVA3158	BLACK	4-0	45	3/8 CIRCLE		13	MICRO POINT	12
195	ADVA3159	BLACK	4-0	45	3/8 CIRCLE		16	MICRO POINT	12
196	ADVA3160	BLACK	4-0	45	3/8 CIRCLE		19	MICRO POINT	12
197	ADVA3161	BLACK	4-0	75	3/8 CIRCLE		19	MICRO POINT	12
198	ADVA3318	BLACK	4-0	70	3/8 CIRCLE		16	-	12
199	ADVA3319	BLACK	4-0	90	3/8 CIRCLE		16	PRIME	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
200	ADVA3319S	BLACK	4-0	35	3/8 CIRCLE		16	PRIME	12
201	ADVA3326	BLACK	4-0	38	3/8 CIRCLE		10	-	12
202	ADVA3356	BLACK	4-0	45	3/8 CIRCLE		19	PRIME	12
203	ADVA3357	BLACK	4-0	60	3/8 CIRCLE		19	PRIME	12
204	ADVA3567	BLACK	4-0	70	3/8 CIRCLE		19	PRIME	12
205	ADVA7419	BLACK	4-0	75	3/8 CIRCLE		26	-	12
206	ADVA1618	BLACK	5-0	45	3/8 CIRCLE		19	PRIME	12
207	ADVA1618T	BLACK	5-0	45	3/8 CIRCLE		19	PRIME	24
208	ADVA3003	BLACK	5-0	90	3/8 CIRCLE		19	PRIME	12
209	ADVA3016	BLACK	5-0	45	3/8 CIRCLE		16	-	12
210	ADVA3018	BLACK	5-0	45	1/2 CIRCLE		15	SLIM BLADE	12
211	ADVA3047	BLACK	5-0	45	3/8 CIRCLE		11	-	12
212	ADVA3048	BLACK	5-0	75	3/8 CIRCLE		16		12
213	ADVA3052	BLACK	5-0	75	3/8 CIRCLE		16	-	12
214	ADVA3060	BLACK	5-0	90	3/8 CIRCLE		16	-	12
215	ADVA3081	BLACK	5-0	76	3/8 CIRCLE		15	-	12
216	ADVA3092	BLACK	5-0	45	3/8 CIRCLE		12	-	12
217	ADVA3102	BLACK	5-0	75	3/8 CIRCLE		19	PRIME	12
218	ADVA3103	BLACK	5-0	75	3/8 CIRCLE		19	PRIME	12
219	ADVA3162	BLACK	5-0	45	3/8 CIRCLE		11	MICRO POINT	12
220	ADVA3163	BLACK	5-0	45	3/8 CIRCLE		13	MICRO POINT	12
221	ADVA3164	BLACK	5-0	45	3/8 CIRCLE		16	MICRO POINT	12
222	ADVA3165	BLACK	5-0	45	3/8 CIRCLE		19	MICRO POINT	12
223	ADVA3315	BLACK	5-0	35	3/8 CIRCLE		16	PRIME	12
224	ADVA3316	BLACK	5-0	90	3/8 CIRCLE		12	PRIME	12
225	ADVA3317	BLACK	5-0	70	3/8 CIRCLE		12	-	12
226	ADVA3353	BLACK	5-0	60	3/8 CIRCLE		13	PRIME	12
227	ADVA3667	BLACK	5-0	75	3/8 CIRCLE		11	SLIM BLADE	12
228	ADVA3726	BLACK	5-0	70	3/8 CIRCLE		15	-	12
229	ADVA3752	BLACK	5-0	70	3/8 CIRCLE		13	PRIME	12
230	ADVA1615T	BLACK	6-0	45	3/8 CIRCLE		16	PRIME	24
231	ADVA3019	BLACK	6-0	35	1/2 CIRCLE		16	-	12
232	ADVA3049	BLACK	6-0	45	3/8 CIRCLE		12	-	12
233	ADVA3093	BLACK	6-0	45	3/8 CIRCLE		15	-	12
234	ADVA3115	BLACK	6-0	75	3/8 CIRCLE		9	PRIME	12
235	ADVA3166	BLACK	6-0	45	3/8 CIRCLE		11	MICRO POINT	12
236	ADVA3167	BLACK	6-0	45	3/8 CIRCLE		13	MICRO POINT	12
237	ADVA3182	BLACK	6-0	45	3/8 CIRCLE		16	-	12
238	ADVA3320	BLACK	6-0	70	3/8 CIRCLE		10	-	12
239	ADVA3331	BLACK	6-0	70	3/8 CIRCLE		26	-	12
240	ADVA3351	BLACK	6-0	70	3/8 CIRCLE		13	PRIME	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
241	ADVA3645	BLACK	6-0	45	3/8 CIRCLE		12	-	12
242	ADVA3668	BLACK	6-0	45	3/8 CIRCLE		16	-	12
243	ADVA3668H	BLACK	6-0	45	3/8 CIRCLE		16	-	36
244	ADVA3945	BLACK	6-0	45	3/8 CIRCLE		9	-	12
245	ADVAW511SA	BLACK	6-0	60	3/8 CIRCLE		13	-	12
246	ADVA2712	BLACK	7-0	60	3/8 CIRCLE		6	-	12
247	ADVA3045	BLACK	7-0	30	3/8 CIRCLE		8	-	12
248	ADVA3109	BLACK	7-0	45	3/8 CIRCLE		11	-	12
249	ADVA2812	BLACK	8-0	13	3/8 CIRCLE		5	-	12
250	ADVA2808G	BLACK	8-0	13	3/8 CIRCLE		6.5	MICRO POINT	12
251	ADVA2908	BLACK	8-0	13	3/8 CIRCLE		6.5	ADVABLACK	12
252	ADVA3020	BLACK	8-0	45	1/2 CIRCLE	 P	8	MICRO POINT	12
253	ADVA3025	BLACK	8-0	13	1/2 CIRCLE		6.5	-	12
254	ADVA3113	BLACK	8-0	30	1/2 CIRCLE		7	-	12
255	ADVA3322	BLACK	8-0	38	3/8 CIRCLE		6	-	12
256	ADVA3343	BLACK	8-0	38	3/8 CIRCLE		8	MICRO POINT	12
257	ADVA3708	BLACK	8-0	38	3/8 CIRCLE		6	MICRO POINT	12
258	ADVA3708SP	BLACK	8-0	38	3/8 CIRCLE		6	-	12
259	ADVA3860	BLACK	8-0	60	3/8 CIRCLE		9	MICROPOINT	12
260	ADVA2829G	BLACK	9-0	13	3/8 CIRCLE		5	MICRO POINT	12
261	ADVA2829L	BLACK	9-0	70	3/8 CIRCLE		5	MICRO POINT	12
262	ADVA2871	BLACK	9-0	13	3/8 CIRCLE		4	-	12
263	ADVA3023	BLACK	9-0	13	1/2 CIRCLE		5	-	12
264	ADVA3024	BLACK	9-0	13	1/2 CIRCLE		4	MICRO POINT	12
265	ADVA3181	BLACK	9-0	15	3/8 CIRCLE		6	-	12
266	ADVA3709	BLACK	9-0	38	3/8 CIRCLE		6	MICRO POINT	12
267	ADVA3715	BLACK	9-0	38	3/8 CIRCLE		6	-	12
268	ADVA3913	BLACK	9-0	13	1/2 CIRCLE		6.5	ADVABLACK	12
269	ADVA2850G	BLACK	10-0	13	3/8 CIRCLE		5	MICRO POINT	12
270	ADVA2870G	BLACK	10-0	13	3/8 CIRCLE		4	MICRO POINT	12
271	ADVA3026	BLACK	10-0	13	1/2 CIRCLE		5	-	12
272	ADVA3031	BLACK	10-0	13	1/2 CIRCLE		4	MICRO POINT	12
273	ADVA3070	BLACK	10-0	15	3/8 CIRCLE		5	MICRO POINT	12
274	ADVA3111	BLACK	10-0	38	3/8 CIRCLE		8	-	36
275	ADVA3314	BLACK	10-0	25	3/8 CIRCLE		4	MICRO POINT	12
276	ADVA3718	BLACK	10-0	38	3/8 CIRCLE		6	-	12
277	ADVA3720	BLACK	10-0	38	3/8 CIRCLE		6	-	12
278	ADVA2729	BLACK	11-0	25	3/8 CIRCLE		6	MICRO POINT	12

 Round Bodied

 Reverse Cutting

 Cutting

 Straight Cutting

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
279	ADVA689	BLACK	1	75	3/8 CIRCLE			90	12
280	ADVA3008	BLACK	1	75	3/8 CIRCLE			90	12
281	ADVA3040	BLACK	2-0	75	1/2 CIRCLE			30	12
282	ADVA8170	BLACK	8-0	45	3/8 CIRCLE			9	12
283	ADVA3324	BLACK	8-0	38	3/8 CIRCLE			6	12
284	ADVA2800G	BLACK	9-0	25	3/8 CIRCLE			6.5	12
285	ADVA3046	BLACK	9-0	30	3/8 CIRCLE			6	12
286	ADVA3067	BLACK	9-0	45	3/8 CIRCLE			7	12
287	ADVA3714	BLACK	9-0	38	3/8 CIRCLE			6	12
288	ADVA3022	BLACK	10-0	30	1/2 CIRCLE			6	12
289	ADVA3058	BLACK	10-0	30	3/8 CIRCLE			6	12
290	ADVA3310	BLACK	10-0	30	3/8 CIRCLE			6.5	12
291	ADVA3311	BLACK	10-0	25	3/8 CIRCLE			4	12
292	ADVA3313	BLACK	10-0	5	3/8 CIRCLE			4	12
293	ADVA3066	BLACK	10-0	45	3/8 CIRCLE			6	12
294	ADVA3717	BLACK	10-0	20	3/8 CIRCLE			6	12
295	ADVA3719	BLACK	10-0	38	3/8 CIRCLE			6	12
296	ADVA2741	BLACK	11-0	25	3/8 CIRCLE			6	12

 Round Bodied

 Reverse Cutting

 Cutting

 Spatulated Micro Point

Product Description


Composition	Polyamide
Coating	None
Sizes	U.S.P. Size 1 to 11-0
Colour	Black
Absorption Profile	Non-Absorbable
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Trocar, Reverse Cutting, Taper Cut, Blunt Point, Oval Round Bodied, Diamond Point & Spatula Micro Point needles

Tensile Strength

Permanent, however progressive hydrolysis of the suture in-vivo may result in gradual loss of its tensile strength over time.

115% more knot pull tensile strength vs. U.S.P. limit ensures high tensile strength for secure retention suturing


Knot pull tensile strength (KgF)¹




REFERENCE _____

1 : Data on file with AMS - Bench test done on Nylon, U.S.P. Size - 2-0, 45 mm Reverse cutting needle

Monofilament Material



Premium Needles for consistent Superior Performance



Minimal tissue trauma enables faster recovery



P-21-22, 25-26 & 34-35, GIDC, Manjusar, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India
T : (91) 2667 671302 | amsltd.com | info@amsltd.com

Singapore Office : 111, Somerset Road, 07-08, 111 Somerset, Singapore - 238164.
T: (65) 636910 80 / 89 | E: info@amsltd.com | amsltd.com