



ADVALENE Polypropylene

Monofilament Synthetic Non-Absorbable Surgical Suture

**Reliable performance for confidence in
contaminated wound closure**



**ADVA TRAY™ - the confidence to
suture with ease**

With exceptional strength comes confidence that lasts

Product Description

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVA805	BLUE	0	70	1/2 CIRCLE		25	-	12
2	ADVA823	BLUE	0	45	3/8 CIRCLE		26	PRIME	12
3	ADVAP824	BLUE	0	100	1/2 CIRCLE		36	-	12
4	ADVA830	BLUE	0	70	1/2 CIRCLE		30	HEAVY	12
5	ADVA833	BLUE	0	150	1/2 CIRCLE		40	-	12
6	ADVAP833	BLUE	0	90	1/2 CIRCLE		50	-	12
7	ADVA842	BLUE	0	70	1/2 CIRCLE		30	-	12
8	ADVA846	BLUE	0	90	1/2 CIRCLE		30	-	12
9	ADVA846SP	BLUE	0	75	1/2 CIRCLE		37	-	12
10	ADVA894	BLUE	0	90	1/2 CIRCLE		30	HEAVY	12
11	ADVA8014	BLUE	0	75	1/2 CIRCLE		26	-	12
12	ADVA8040	BLUE	0	70	1/2 CIRCLE		36	-	12
13	ADVA8041	BLUE	0	70	1/2 CIRCLE		40	-	12
14	ADVA8074	BLUE	0	100	1/2 CIRCLE		40	HEAVY	12
15	ADVA8075	BLUE	0	100	1/2 CIRCLE		40	-	12
16	ADVA8162	BLUE	0	100	1/2 CIRCLE		40	HEAVY	12
17	ADVA8165	BLUE	0	100	J TYPE		31	-	12
18	ADVA8166	BLUE	0	100	1/2 CIRCLE		48	-	12
19	ADVA8174	BLUE	0	75	1/2 CIRCLE		30	-	12
20	ADVA8226	BLUE	0	70	3/8 CIRCLE		12	-	12
21	ADVA8235	BLUE	0	100	1/2 CIRCLE		30	-	12
22	ADVA8424	BLUE	0	75	1/2 CIRCLE		36	-	12
23	ADVA8690	BLUE	0	76	3/8 CIRCLE		36	-	12
24	ADVA831	BLUE	1	100	1/2 CIRCLE		60	-	12
25	ADVA834	BLUE	1	150	1/2 CIRCLE		48	HEAVY	12
26	ADVA834H	BLUE	1	150	1/2 CIRCLE		40	HEAVY	12
27	ADVAP834	BLUE	1	100	1/2 CIRCLE		50	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

 LP - Loop

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
28	ADVA840	BLUE	1	70	1/2 CIRCLE		30	HEAVY	12
29	ADVA843	BLUE	1	70	1/2 CIRCLE		40	HEAVY	12
30	ADVA883	BLUE	1	100	1/2 CIRCLE		45	HEAVY	12
31	ADVA896	BLUE	1	100	1/2 CIRCLE		30	HEAVY	12
32	ADVA8035	BLUE	1	100	1/2 CIRCLE		40	-	12
33	ADVA8038	BLUE	1	90	1/2 CIRCLE		75	-	12
34	ADVA8039	BLUE	1	70	1/2 CIRCLE		40	-	12
35	ADVA8044	BLUE	1	75	3/8 CIRCLE		40	-	12
36	ADVA8045	BLUE	1	90	1/2 CIRCLE		45	-	12
37	ADVA8046	BLUE	1	100	1/2 CIRCLE		50	-	12
38	ADVA8054	BLUE	1	75	1/2 CIRCLE		37	-	12
39	ADVA8073	BLUE	1	150	1/2 CIRCLE	 LP	40	-	12
40	ADVA8083	BLUE	1	150	1/2 CIRCLE		48	HEAVY	12
41	ADVA8084	BLUE	1	75	3/8 CIRCLE		36	-	12
42	ADVA8091	BLUE	1	90	1/2 CIRCLE		45	-	12
43	ADVA8116	BLUE	1	90	1/2 CIRCLE		24	-	12
44	ADVA8118	BLUE	1	90	1/2 CIRCLE		30	-	12
45	ADVA8126	BLUE	1	70	3/8 CIRCLE		12	-	12
46	ADVA8129	BLUE	1	75	1/2 CIRCLE		40	-	12
47	ADVA8173	BLUE	1	100	1/2 CIRCLE		45	-	12
48	ADVA8178	BLUE	1	100	1/2 CIRCLE		45	-	12
49	ADVA3846	BLUE	2	3 X 50	3/8 CIRCLE		75	HEAVY	12
50	ADVA8090	BLUE	2	76	1/2 CIRCLE		50	-	12
51	ADVA8130	BLUE	2	45	1/2 CIRCLE		26	-	12
52	ADVA806	BLUE	2-0	90	1/2 CIRCLE		26	-	12
53	ADVA807	BLUE	2-0	70	1/2 CIRCLE		26	-	12
54	ADVA819	BLUE	2-0	45	3/8 CIRCLE		15	-	12
55	ADVA833SL	BLUE	2-0	90	1/2 CIRCLE		30	-	12
56	ADVA841	BLUE	2-0	70	1/2 CIRCLE		30	-	12
57	ADVA841SP	BLUE	2-0	90	1/2 CIRCLE		20	-	12
58	ADVA844	BLUE	2-0	90	1/2 CIRCLE		30	-	12
59	ADVA850	BLUE	2-0	90	1/2 CIRCLE		36	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
60	ADVA863	BLUE	2-0	70	3/8 CIRCLE		45	-	12
61	ADVA8005	BLUE	2-0	90	1/2 CIRCLE		26	-	12
62	ADVA8008	BLUE	2-0	70	3/8 CIRCLE		30	-	12
63	ADVA8008S	BLUE	2-0	45	3/8 CIRCLE		30	-	12
64	ADVA8009	BLUE	2-0	75	1/2 CIRCLE		30	-	12
65	ADVA8026L	BLUE	2-0	60	3/8 CIRCLE		26	PRIME	12
66	ADVA8026T	BLUE	2-0	45	3/8 CIRCLE		26	PRIME	24
67	ADVA8051	BLUE	2-0	45	3/8 CIRCLE		24	-	12
68	ADVA8076	BLUE	2-0	75	3/8 CIRCLE		20	-	12
69	ADVA8101	BLUE	2-0	75	3/8 CIRCLE		30	-	12
70	ADVA8111	BLUE	2-0	45	3/8 CIRCLE		26	PRIME	12
71	ADVA8115	BLUE	2-0	90	1/2 CIRCLE		24	-	12
72	ADVA8119	BLUE	2-0	90	1/2 CIRCLE		25	-	12
73	ADVA8121	BLUE	2-0	90	1/2 CIRCLE		30	-	12
74	ADVA8153	BLUE	2-0	75	1/2 CIRCLE		30	-	12
75	ADVA8167	BLUE	2-0	100	3/8 CIRCLE		40	-	12
76	ADVA8168	BLUE	2-0	90	STRAIGHT		60	-	12
77	ADVA8171	BLUE	2-0	75	3/8 CIRCLE		24	-	12
78	ADVA8245	BLUE	2-0	45	1/2 CIRCLE		25	-	12
79	ADVA8245L	BLUE	2-0	100	1/2 CIRCLE		25	-	12
80	ADVA8251	BLUE	2-0	45	3/8 CIRCLE		13	-	12
81	ADVA8411	BLUE	2-0	75	1/2 CIRCLE		26	-	12
82	ADVA8433	BLUE	2-0	75	1/2 CIRCLE		40	-	12
83	ADVA8571	BLUE	2-0	90	STRAIGHT		55	-	12
84	ADVA8631	BLUE	2-0	75	STRAIGHT		60	-	12
85	ADVA8685	BLUE	2-0	75	3/8 CIRCLE		24	-	12
86	ADVA8689	BLUE	2-0	76	3/8 CIRCLE		24	-	12
87	ADVA8833	BLUE	2-0	76	1/2 CIRCLE		26	-	12
88	ADVA8937S	BLUE	2-0	75	3/8 CIRCLE		27	-	12
89	ADVA018	BLUE	3-0	70	3/8 CIRCLE		22	-	12
90	ADVA018S	BLUE	3-0	45	3/8 CIRCLE		22	-	12
91	ADVA019	BLUE	3-0	70	3/8 CIRCLE		16	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
92	ADVAP422	BLUE	3-0	60	3/8 CIRCLE		26	-	12
93	ADVA0621	BLUE	3-0	90	STRAIGHT		60	-	12
94	ADVA800	BLUE	3-0	90	3/8 CIRCLE		25	PRIME	12
95	ADVA802	BLUE	3-0	70	3/8 CIRCLE		25	-	12
96	ADVA802SP	BLUE	3-0	75	3/8 CIRCLE		13	-	12
97	ADVA820	BLUE	3-0	45	3/8 CIRCLE		25	PRIME	12
98	ADVA821	BLUE	3-0	60	3/8 CIRCLE		25	-	12
99	ADVA822L	BLUE	3-0	75	1/2 CIRCLE		15	-	12
100	ADVA825	BLUE	3-0	90	1/2 CIRCLE		25	-	12
101	ADVA838	BLUE	3-0	70	1/2 CIRCLE		25	-	12
102	ADVA860	BLUE	3-0	75	STRAIGHT		60	-	12
103	ADVA887	BLUE	3-0	70	1/2 CIRCLE		20	-	12
104	ADVA8001	BLUE	3-0	70	1/2 CIRCLE		22	-	12
105	ADVA8003	BLUE	3-0	70	1/2 CIRCLE		36	-	12
106	ADVA8004	BLUE	3-0	70	1/2 CIRCLE		40	-	12
107	ADVA8006	BLUE	3-0	45	3/8 CIRCLE		26	PRIME	12
108	ADVA8021	BLUE	3-0	75	1/2 CIRCLE		22	-	12
109	ADVA8021T	BLUE	3-0	45	3/8 CIRCLE		26	P Curved	24
110	ADVA8042	BLUE	3-0	70	1/2 CIRCLE		24	-	12
111	ADVA8052	BLUE	3-0	75	3/8 CIRCLE		24	-	12
112	ADVA8061	BLUE	3-0	75	1/2 CIRCLE		16	-	12
113	ADVA8077	BLUE	3-0	75	3/8 CIRCLE		24	-	12
114	ADVA8078	BLUE	3-0	75	1/2 CIRCLE		18	-	12
115	ADVA8110	BLUE	3-0	45	3/8 CIRCLE		26	P Curved	12
116	ADVA8132	BLUE	3-0	70	3/8 CIRCLE		12	-	12
117	ADVA8139	BLUE	3-0	45	3/8 CIRCLE		19	-	12
118	ADVA8140	BLUE	3-0	45	3/8 CIRCLE		24	-	12
119	ADVA8163	BLUE	3-0	90	1/2 CIRCLE		17	-	12
120	ADVA8164	BLUE	3-0	45	3/8 CIRCLE		19	PRIME	12
121	ADVA8172	BLUE	3-0	75	3/8 CIRCLE		24	-	12
122	ADVA8175	BLUE	3-0	75	3/8 CIRCLE		22	-	12
123	ADVA8558SL	BLUE	3-0	75	1/2 CIRCLE		17	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
124	ADVA8675	BLUE	3-0	45	3/8 CIRCLE		30	-	12
125	ADVA8684	BLUE	3-0	76	3/8 CIRCLE		24	-	12
126	ADVA8770	BLUE	3-0	60	1/2 CIRCLE		22	-	12
127	ADVA8771XL	BLUE	3-0	90	1/2 CIRCLE		22	-	12
128	ADVA8832	BLUE	3-0	75	1/2 CIRCLE		26	LONG TAPER	12
129	ADVA804	BLUE	4-0	70	3/8 CIRCLE		15	-	12
130	ADVA815	BLUE	4-0	45	3/8 CIRCLE		16	-	12
131	ADVA816	BLUE	4-0	60	3/8 CIRCLE		16	PRIME	12
132	ADVA848	BLUE	4-0	90	1/2 CIRCLE		16	LONG TAPER	12
133	ADVA849	BLUE	4-0	70	1/2 CIRCLE		16	CV	12
134	ADVA851	BLUE	4-0	70	3/8 CIRCLE		16	PRIME	12
135	ADVA857	BLUE	4-0	90	1/2 CIRCLE		17	-	12
136	ADVA870	BLUE	4-0	70	3/8 CIRCLE		16	PRIME	12
137	ADVA872	BLUE	4-0	90	3/8 CIRCLE		25	-	12
138	ADVA879	BLUE	4-0	45	3/8 CIRCLE		12	-	12
139	ADVA879L	BLUE	4-0	75	3/8 CIRCLE		12	-	12
140	ADVA899	BLUE	4-0	70	1/2 CIRCLE		17	-	12
141	ADVA3986	BLUE	4-0	75	1/2 CIRCLE		22	-	12
142	ADVA8007	BLUE	4-0	45	3/8 CIRCLE		19	PRIME	12
143	ADVA8043	BLUE	4-0	75	1/2 CIRCLE		26	-	12
144	ADVA8047	BLUE	4-0	60	3/8 CIRCLE		16	PRIME	12
145	ADVA8053	BLUE	4-0	75	3/8 CIRCLE		19	-	12
146	ADVA8058	BLUE	4-0	90	1/2 CIRCLE		13	-	12
147	ADVA8064	BLUE	4-0	45	3/8 CIRCLE		25	PRIME	12
148	ADVA8079	BLUE	4-0	150	3/8 CIRCLE		19	-	12
149	ADVA8092	BLUE	4-0	90	3/8 CIRCLE		26	-	12
150	ADVA8103	BLUE	4-0	60	3/8 CIRCLE		10	-	12
151	ADVA8112	BLUE	4-0	90	1/2 CIRCLE		19	-	12
152	ADVA8120	BLUE	4-0	75	1/2 CIRCLE		30	-	12
153	ADVA8122	BLUE	4-0	90	1/2 CIRCLE		30	-	12
154	ADVA8127	BLUE	4-0	70	3/8 CIRCLE		12	-	12
155	ADVA8137	BLUE	4-0	45	3/8 CIRCLE		12	-	12

 Round Bodied Reverse Cutting Cutting Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
156	ADVA8143	BLUE	4-0	45	3/8 CIRCLE		19	-	12
157	ADVA8149	BLUE	4-0	60	3/8 CIRCLE		17	-	12
158	ADVA8157	BLUE	4-0	90	1/2 CIRCLE		26	LONG TAPER	12
159	ADVA8169	BLUE	4-0	45	3/8 CIRCLE		16	SLIM BLADE	12
160	ADVA8682	BLUE	4-0	45	3/8 CIRCLE		19	PRIME	12
161	ADVA8871	BLUE	4-0	75	1/2 CIRCLE		17	-	12
162	ADVA826	BLUE	5-0	70	3/8 CIRCLE		15	LONG TAPER	12
163	ADVA826S	BLUE	5-0	45	3/8 CIRCLE		15	LONG TAPER	12
164	ADVA862	BLUE	5-0	70	3/8 CIRCLE		13	PRIME	12
165	ADVA880	BLUE	5-0	70	3/8 CIRCLE		12	PRIME	12
166	ADVA881	BLUE	5-0	70	3/8 CIRCLE		16	CV	12
167	ADVA881SP	BLUE	5-0	75	3/8 CIRCLE		18	-	12
168	ADVA889	BLUE	5-0	90	3/8 CIRCLE		12	PRIME	12
169	ADVA8010	BLUE	5-0	75	1/2 CIRCLE		18	-	12
170	ADVA8031	BLUE	5-0	75	1/2 CIRCLE		13	320 MICRONS	12
171	ADVA8048	BLUE	5-0	60	3/8 CIRCLE		16	PRIME	12
172	ADVA8055	BLUE	5-0	90	1/2 CIRCLE		17	PRIME	12
173	ADVA8065	BLUE	5-0	45	3/8 CIRCLE		20	PRIME	12
174	ADVA8080	BLUE	5-0	75	3/8 CIRCLE		12	-	12
175	ADVA8113	BLUE	5-0	90	1/2 CIRCLE		19	-	12
176	ADVA8123	BLUE	5-0	90	1/2 CIRCLE		30	-	12
177	ADVA8138	BLUE	5-0	45	3/8 CIRCLE		19	-	12
178	ADVA8144	BLUE	5-0	90	1/2 CIRCLE		17	-	12
179	ADVA8145	BLUE	5-0	45	3/8 CIRCLE		13	PRIME	12
180	ADVA8150	BLUE	5-0	60	3/8 CIRCLE		17	-	12
181	ADVA8159	BLUE	5-0	70	3/8 CIRCLE		15	SLIM BLADE	12
182	ADVA8170	BLUE	5-0	45	3/8 CIRCLE		15	PRIME	12
183	ADVA8180XL	BLUE	5-0	75	3/8 CIRCLE		19	-	12
184	ADVA8527	BLUE	5-0	45	3/8 CIRCLE		16	PRIME	12
185	ADVA8635	BLUE	5-0	45	3/8 CIRCLE		16	PRIME	12
186	ADVA8872	BLUE	5-0	45	3/8 CIRCLE		13	PRIME	12
187	ADVA8872T	BLUE	5-0	45	3/8 CIRCLE		13	PRIME	24

 Round Bodied Reverse Cutting Cutting Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
188	ADVA8882	BLUE	5-0	45	3/8 CIRCLE		16	PRIME	12
189	ADVA8882T	BLUE	5-0	45	3/8 CIRCLE		16	PRIME	24
190	ADVA861	BLUE	6-0	75	3/8 CIRCLE		8	-	12
191	ADVA868	BLUE	6-0	75	3/8 CIRCLE		12	-	12
192	ADVA868SP	BLUE	6-0	60	3/8 CIRCLE		8	220 MICRONS	12
193	ADVA878	BLUE	6-0	45	3/8 CIRCLE		12	-	12
194	ADVA878SP	BLUE	6-0	45	3/8 CIRCLE		15	-	12
195	ADVAP875	BLUE	6-0	75	1/2 CIRCLE		8	-	12
196	ADVA8060	BLUE	6-0	60	3/8 CIRCLE		10	-	12
197	ADVA8066	BLUE	6-0	75	3/8 CIRCLE		11	PRIME	12
198	ADVA8098	BLUE	6-0	60	3/8 CIRCLE		13	-	12
199	ADVA8114	BLUE	6-0	90	1/2 CIRCLE		19	-	12
200	ADVA8117	BLUE	6-0	45	3/8 CIRCLE		16	PRIME	12
201	ADVA8141	BLUE	6-0	45	3/8 CIRCLE		13	-	12
202	ADVA8142	BLUE	6-0	60	3/8 CIRCLE		11	-	12
203	ADVA8158	BLUE	6-0	90	3/8 CIRCLE		13	ADVABLACK	12
204	ADVA8176	BLUE	6-0	75	3/8 CIRCLE		12	-	12
205	ADVA8221	BLUE	6-0	45	3/8 CIRCLE		13	PRIME	12
206	ADVA8680	BLUE	6-0	45	3/8 CIRCLE		16	-	12
207	ADVA8868T	BLUE	6-0	45	3/8 CIRCLE		11	PRIME	24
208	ADVA811SA	BLUE	7-0	60	3/8 CIRCLE		9.5	-	12
209	ADVAP837	BLUE	7-0	75	1/2 CIRCLE		13	-	12
210	ADVA2076SA	BLUE	7-0	75	3/8 CIRCLE		6	-	12
211	ADVA2078L	BLUE	7-0	75	1/2 CIRCLE		8	-	12
212	ADVA8068	BLUE	7-0	60	3/8 CIRCLE		8	-	12
213	ADVA8068H	BLUE	7-0	60	3/8 CIRCLE		8	-	36
214	ADVA8069	BLUE	7-0	75	3/8 CIRCLE		8	-	12
215	ADVA8088SP	BLUE	7-0	75	3/8 CIRCLE		8	-	12
216	ADVA8107	BLUE	7-0	90	1/2 CIRCLE		13	-	12
217	ADVA900SP	BLUE	8-0	38	3/8 CIRCLE		6	-	12
218	ADVA8033	BLUE	8-0	13	3/8 CIRCLE		5	-	12
219	ADVA8070	BLUE	8-0	60	3/8 CIRCLE		8	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
220	ADVA8102	BLUE	8-0	60	1/4 CIRCLE		6	-	12
221	ADVA8131	BLUE	8-0	45	1/2 CIRCLE		8	-	12
222	ADVA2680	BLUE	9-0	13	3/8 CIRCLE		4.7	-	12
223	ADVA2780	BLUE	9-0	13	3/8 CIRCLE		5	-	12
224	ADVA3709SP	BLUE	9-0	38	3/8 CIRCLE		6	MICRO POINT	12
225	ADVA8022	BLUE	9-0	23	3/8 CIRCLE		6.5	-	12
226	ADVA1613	BLUE	10-0	13	3/8 CIRCLE		3.8	-	12
227	ADVA2790	BLUE	10-0	13	3/8 CIRCLE		4	-	12
228	ADVA2794	BLUE	10-0	13	3/8 CIRCLE		5	-	12
229	ADVA8049	BLUE	10-0	25	3/8 CIRCLE		6	-	12
230	ADVA8154	BLUE	10-0	23	1/2 CIRCLE		16	MICRO POINT	12

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
231	ADVA824	BLUE	0	90	1/2 CIRCLE	 LT	 LT	26	12
232	ADVA8056	BLUE	0	75	1/2 CIRCLE			13	12
233	ADVA8179	BLUE	0	100	1/2 CIRCLE			36	12
234	ADVA808	BLUE	2-0	70	1/2 CIRCLE			25	12
235	ADVA832	BLUE	2-0	70	STRAIGHT			30	12
236	ADVA837	BLUE	2-0	90	1/2 CIRCLE	 CV	 CV	17	12
237	ADVA877	BLUE	2-0	90	1/2 CIRCLE			25	12
238	ADVA898	BLUE	2-0	70	1/2 CIRCLE			17	12
239	ADVA8128	BLUE	2-0	70	1/2 CIRCLE			20	12
240	ADVA8526	BLUE	2-0	90	1/2 CIRCLE	 LT	 LT	31	12
241	ADVA8526SL	BLUE	2-0	75	1/2 CIRCLE			22	12
242	ADVA8876	BLUE	2-0	90	1/2 CIRCLE			24	12
243	ADVA8886	BLUE	2-0	90	1/2 CIRCLE			26	12
244	ADVA8922	BLUE	2-0	90	3/8 CIRCLE			22	12
245	ADVA8937	BLUE	2-0	90	1/2 CIRCLE	 CV	 CV	17	12
246	ADVA8977	BLUE	2-0	90	1/2 CIRCLE			26	12
247	ADVA818	BLUE	3-0	60	3/8 CIRCLE	 LT	 LT	13	12
248	ADVA822	BLUE	3-0	70	1/2 CIRCLE	 LT	 LT	25	12
249	ADVA836	BLUE	3-0	90	1/2 CIRCLE			17	12
250	ADVA852	BLUE	3-0	90	1/2 CIRCLE	 LT	 LT	25	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
251	ADVA897	BLUE	3-0	70	1/2 CIRCLE			17	12
252	ADVA8002	BLUE	3-0	70	3/8 CIRCLE			17	12
253	ADVA8002S	BLUE	3-0	75	3/8 CIRCLE			17	12
254	ADVA8018	BLUE	3-0	90	3/8 CIRCLE			31	12
255	ADVA8036	BLUE	3-0	90	1/2 CIRCLE			25	12
256	ADVA8057	BLUE	3-0	75	1/2 CIRCLE			13	12
257	ADVA8085	BLUE	3-0	75	1/2 CIRCLE			18	12
258	ADVA8100	BLUE	3-0	90	1/2 CIRCLE			26	12
259	ADVA8108	BLUE	3-0	90	1/2 CIRCLE			25	12
260	ADVA8124	BLUE	3-0	70	1/2 CIRCLE			25	12
261	ADVA8146	BLUE	3-0	75	1/2 CIRCLE			20	12
262	ADVA8332	BLUE	3-0	90	1/2 CIRCLE			30	12
263	ADVA8422	BLUE	3-0	90	1/2 CIRCLE			24	12
264	ADVA8522	BLUE	3-0	90	1/2 CIRCLE			26	12
265	ADVA8525	BLUE	3-0	90	1/2 CIRCLE	 LT	 LT	31	12
266	ADVA8558	BLUE	3-0	90	1/2 CIRCLE			17	12
267	ADVA8667	BLUE	3-0	90	1/2 CIRCLE			16	12
268	ADVA8771	BLUE	3-0	75	1/2 CIRCLE			22	12
269	ADVA8936	BLUE	3-0	90	1/2 CIRCLE	 CV	 CV	17	12
270	ADVA8976	BLUE	3-0	90	1/2 CIRCLE			26	12
271	ADVA835	BLUE	4-0	70	1/2 CIRCLE	 CV	 CV	17	12
272	ADVA835SL	BLUE	4-0	45	1/2 CIRCLE			17	12
273	ADVA845	BLUE	4-0	90	1/2 CIRCLE	 CV	 CV	25	12
274	ADVA853	BLUE	4-0	90	1/2 CIRCLE			17	12
275	ADVA867	BLUE	4-0	90	1/2 CIRCLE			25	12
276	ADVA870DA	BLUE	4-0	70	3/8 CIRCLE	 P	 P	16	12
277	ADVA7535PL3	BLUE	4-0	2 X 75	1/2 CIRCLE		 CV	22	12
278	ADVA7588PL3	BLUE	4-0	2 X 90	1/2 CIRCLE		 CV	26	12
279	ADVA8015	BLUE	4-0	90	1/2 CIRCLE			22	12
280	ADVA8019	BLUE	4-0	90	1/2 CIRCLE			40	12
281	ADVA8020	BLUE	4-0	90	1/2 CIRCLE			36	12
282	ADVA8023	BLUE	4-0	90	1/2 CIRCLE			20	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
283	ADVA8026	BLUE	4-0	90	1/2 CIRCLE			16	12
284	ADVA8034	BLUE	4-0	90	1/2 CIRCLE			16	12
285	ADVA8086	BLUE	4-0	60	1/2 CIRCLE			13	12
286	ADVA8094	BLUE	4-0	90	1/2 & 3/8 CIRCLE	 ADVABLACK	 ADVABLACK	17	12
287	ADVA8099	BLUE	4-0	75	3/8 CIRCLE			8	12
288	ADVA8125	BLUE	4-0	70	1/2 CIRCLE			25	12
289	ADVA8134	BLUE	4-0	90	3/8 CIRCLE			20	12
290	ADVA8160	BLUE	4-0	90	1/2 CIRCLE			16	12
291	ADVA8204	BLUE	4-0	60	1/2 CIRCLE	 LT	 LT	13	12
292	ADVA8204PL3	BLUE	4-0	60	1/2 CIRCLE	 CV	 CV	13	12
293	ADVA8354	BLUE	4-0	70	3/8 CIRCLE			16	12
294	ADVA8354L	BLUE	4-0	70	3/8 CIRCLE			16	12
295	ADVA8521	BLUE	4-0	90	1/2 CIRCLE	 LT	 LT	26	12
296	ADVA8521S	BLUE	4-0	70	1/2 CIRCLE			25	12
297	ADVA8521T	BLUE	4-0	90	1/2 CIRCLE	 LT	 LT	26	12
298	ADVA8526CU	BLUE	4-0	75	3/8 CIRCLE			22	12
299	ADVA8526DA	BLUE	4-0	75	1/2 CIRCLE			22	12
300	ADVA8581	BLUE	4-0	90	3/8 CIRCLE			17	12
301	ADVA8871DA	BLUE	4-0	75	1/2 CIRCLE			17	12
302	ADVA8557	BLUE	4-0	90	1/2 CIRCLE	 LT	 LT	17	12
303	ADVA8557PL3	BLUE	4-0	90	1/2 CIRCLE	 CV LONG TAPER	 CV LONG TAPER	17	12
304	ADVA8665	BLUE	4-0	75	1/2 CIRCLE			13	12
305	ADVA8761	BLUE	4-0	90	1/2 CIRCLE	 LT	 LT	20	12
306	ADVA8845	BLUE	4-0	90	1/2 CIRCLE			26	12
307	ADVA8884T	BLUE	4-0	45	1/2 CIRCLE			16 & 26	24
308	ADVA8935	BLUE	4-0	90	1/2 CIRCLE			17	12
309	ADVA8935PL3	BLUE	4-0	90	1/2 CIRCLE	 CV	 CV	17	12
310	ADVA810	BLUE	5-0	75	1/2 CIRCLE	 LT	 LT	13	12
311	ADVA827	BLUE	5-0	70	3/8 CIRCLE	 LT	 LT	12	12
312	ADVA856	BLUE	5-0	90	1/2 CIRCLE	 LT	 LT	16	12
313	ADVA882	BLUE	5-0	70	3/8 CIRCLE	 LT	 LT	16	12
314	ADVA892	BLUE	5-0	70	3/8 CIRCLE	 LT	 LT	20	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
315	ADVA8016	BLUE	5-0	90	3/8 CIRCLE			17	12
316	ADVA8024	BLUE	5-0	90	1/2 CIRCLE			16	12
317	ADVA8027	BLUE	5-0	75	1/2 CIRCLE			13	12
318	ADVA8050	BLUE	5-0	90	1/2 CIRCLE			20	12
319	ADVA8096	BLUE	5-0	90	1/2 CIRCLE			18	12
320	ADVA8133	BLUE	5-0	70	3/8 CIRCLE			12	12
321	ADVA8147	BLUE	5-0	75	1/2 CIRCLE			22	12
322	ADVA8161	BLUE	5-0	90	3/8 CIRCLE			13	12
323	ADVA8205	BLUE	5-0	60	1/2 CIRCLE	 LT	 LT	13	12
324	ADVA8205PL3	BLUE	5-0	60	1/2 CIRCLE	 CV	 CV	13	12
325	ADVA8321	BLUE	5-0	60	3/8 CIRCLE	 ADVABLACK	 ADVABLACK	13	12
326	ADVA8556	BLUE	5-0	90	1/2 CIRCLE			17	12
327	ADVA8556SL	BLUE	5-0	75	1/2 CIRCLE			17	12
328	ADVA8556PL3	BLUE	5-0	90	1/2 CIRCLE	 CV	 CV	17	12
329	ADVA8605PL3	BLUE	5-0	90	1/2 CIRCLE	 CV	 CV	13	12
330	ADVA8664	BLUE	5-0	60	1/2 CIRCLE			13	12
331	ADVA8710	BLUE	5-0	75	1/2 CIRCLE			13	12
332	ADVA8710PL3	BLUE	5-0	75	1/2 CIRCLE	 CV	 CV	13	12
333	ADVA8710S	BLUE	5-0	45	1/2 CIRCLE	 LT	 LT	13	12
334	ADVA8721	BLUE	5-0	75	3/8 CIRCLE			13	12
335	ADVA8721TC	BLUE	5-0	75	3/8 CIRCLE			13	12
336	ADVA8803	BLUE	5-0	60	1/2 CIRCLE			11	12
337	ADVA8803PL3	BLUE	5-0	60	1/2 CIRCLE	 CV	 CV	11	12
338	ADVA8810	BLUE	5-0	60	1/2 CIRCLE			10	12
339	ADVA817	BLUE	6-0	70	3/8 CIRCLE	 P	 P	13	12
340	ADVA868DA	BLUE	6-0	75	3/8 CIRCLE			12	12
341	ADVA8011	BLUE	6-0	75	3/8 CIRCLE			13	12
342	ADVA8025	BLUE	6-0	60	1/2 CIRCLE			13	12
343	ADVA8028	BLUE	6-0	60	1/2 CIRCLE			8	12
344	ADVA8030	BLUE	6-0	60	1/2 CIRCLE			11	12
345	ADVA8032	BLUE	6-0	45	1/2 CIRCLE			8	12
346	ADVA8059	BLUE	6-0	60	3/8 CIRCLE			9.5	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
347	ADVA8062	BLUE	6-0	75	3/8 CIRCLE			10	12
348	ADVA8071	BLUE	6-0	75	3/8 CIRCLE			10	12
349	ADVA8081	BLUE	6-0	75	1/2 CIRCLE			10	12
350	ADVA8087	BLUE	6-0	75	3/8 CIRCLE	 P	 P	13	12
351	ADVA8104	BLUE	6-0	60	3/8 CIRCLE	 ADVAPASS	 ADVAPASS	13	12
352	ADVA8109	BLUE	6-0	60	3/8 CIRCLE			8	12
353	ADVA8135	BLUE	6-0	90	3/8 CIRCLE			13	12
354	ADVA8148	BLUE	6-0	75	1/2 CIRCLE			22	12
355	ADVA8152	BLUE	6-0	45	1/2 CIRCLE			8	12
356	ADVA8177	BLUE	6-0	75	3/8 CIRCLE			21	12
357	ADVA8206PL3	BLUE	6-0	60	1/2 CIRCLE	 CV	 CV	13	12
358	ADVA8207	BLUE	6-0	60	1/2 CIRCLE			10	12
359	ADVA8207PL3	BLUE	6-0	60	1/2 CIRCLE	 CV	 CV	10	12
360	ADVA8216	BLUE	6-0	75	1/2 CIRCLE			13	12
361	ADVA8307	BLUE	6-0	60	3/8 CIRCLE	 ADVABLACK	 ADVABLACK	13	12
362	ADVA8403	BLUE	6-0	60	3/8 CIRCLE			9	12
363	ADVA8597	BLUE	6-0	60	3/8 CIRCLE			11	12
364	ADVA8597SP	BLUE	6-0	60	3/8 CIRCLE			10	12
365	ADVA8681	BLUE	6-0	75	3/8 CIRCLE			13	12
366	ADVA8707	BLUE	6-0	60	3/8 CIRCLE			13	12
367	ADVA8717	BLUE	6-0	90	1/2 CIRCLE			17	12
368	ADVA8718	BLUE	6-0	60	3/8 CIRCLE			13	12
369	ADVA8725T	BLUE	6-0	75	3/8 CIRCLE			13	12
370	ADVA8726	BLUE	6-0	60	3/8 CIRCLE			13	12
371	ADVA8726PL3	BLUE	6-0	60	3/8 CIRCLE	 CV	 CV	13	12
372	ADVA8776	BLUE	6-0	60	3/8 CIRCLE			11	12
373	ADVA8802	BLUE	6-0	60	3/8 CIRCLE			11	12
374	ADVA8804PL3	BLUE	6-0	60	3/8 CIRCLE	 CV	 CV	10	12
375	ADVA8804XLPL3	BLUE	6-0	75	3/8 CIRCLE		 CV	10	12
376	ADVA8807	BLUE	6-0	60	3/8 CIRCLE	 ADVAPASS	 ADVAPASS	9.5	12
377	ADVA8870	BLUE	6-0	45	3/8 CIRCLE	 ADVABLACK	 ADVABLACK	8	12
378	ADVA804DA	BLUE	7-0	60	3/8 CIRCLE			8	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
379	ADVA811	BLUE	7-0	60	3/8 CIRCLE			6	12
380	ADVA811SP	BLUE	7-0	60	3/8 CIRCLE			6.4	12
381	ADVA2013	BLUE	7-0	75	3/8 CIRCLE			13	12
382	ADVA2070L	BLUE	7-0	75	3/8 CIRCLE			10	12
383	ADVA2078	BLUE	7-0	60	1/2 CIRCLE			10	12
384	ADVA2078B	BLUE	7-0	60	COMPOUND CURVE			8	12
385	ADVA2079L	BLUE	7-0	75	3/8 CIRCLE			9	12
386	ADVA8012	BLUE	7-0	60	3/8 CIRCLE			11	12
387	ADVA8017	BLUE	7-0	75	3/8 CIRCLE			9.5	12
388	ADVA8029	BLUE	7-0	60	1/2 CIRCLE			8	12
389	ADVA8037	BLUE	7-0	60	1/2 CIRCLE			8	12
390	ADVA8063	BLUE	7-0	60	3/8 CIRCLE			8	12
391	ADVA8067	BLUE	7-0	60	1/2 CIRCLE			9	12
392	ADVA8078DA	BLUE	7-0	60	3/8 CIRCLE			8	12
393	ADVA8082	BLUE	7-0	75	1/2 CIRCLE			10	12
394	ADVA8088	BLUE	7-0	60	3/8 CIRCLE			9.5	12
395	ADVA8097	BLUE	7-0	75	3/8 CIRCLE			9	12
396	ADVA8151	BLUE	7-0	60	3/8 CIRCLE			13	12
397	ADVA8156	BLUE	7-0	60	3/8 CIRCLE			11	12
398	ADVA8180	BLUE	7-0	75	3/8 CIRCLE			9	12
399	ADVA8188	BLUE	7-0	75	3/8 CIRCLE			7.6	12
400	ADVA8304	BLUE	7-0	60	3/8 CIRCLE			9.5	12
401	ADVA8702	BLUE	7-0	60	3/8 CIRCLE			9	12
402	ADVA8702SP	BLUE	7-0	60	3/8 CIRCLE			9.5	12
403	ADVA8703XL	BLUE	7-0	75	3/8 CIRCLE			9.3	12
404	ADVA8704	BLUE	7-0	60	3/8 CIRCLE			9.5	12
405	ADVA8704PL3	BLUE	7-0	60	3/8 CIRCLE			9.3	12
406	ADVA8725	BLUE	7-0	60	3/8 CIRCLE			13	12
407	ADVA8735	BLUE	7-0	60	3/8 CIRCLE			8	12
408	ADVA8735H	BLUE	7-0	60	3/8 CIRCLE			8	36
409	ADVA8747	BLUE	7-0	60	3/8 CIRCLE			9.5	12
410	ADVA8747H	BLUE	7-0	60	3/8 CIRCLE			9.5	36

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
411	ADVA2088L	BLUE	8-0	75	3/8 CIRCLE			8	12
412	ADVA8013	BLUE	8-0	60	3/8 CIRCLE			6	12
413	ADVA8072	BLUE	8-0	60	3/8 CIRCLE			9.5	12
414	ADVA8089	BLUE	8-0	60	3/8 CIRCLE			8	12
415	ADVA8093	BLUE	8-0	90	3/8 CIRCLE	 ADVABLACK	 ADVABLACK	8	12
416	ADVA8101BL	BLUE	8-0	60	3/8 CIRCLE	 ADVABLACK	 ADVABLACK	8	12
417	ADVA8105	BLUE	8-0	45	3/8 CIRCLE	 ADVAPASS	 ADVAPASS	6.5	12
418	ADVA8106	BLUE	8-0	60	3/8 CIRCLE	 ADVABLACK	 ADVABLACK	8	12
419	ADVA8136	BLUE	8-0	90	1/2 CIRCLE			9	12
420	ADVA8632	BLUE	8-0	75	1/2 CIRCLE			6	12
421	ADVA8703	BLUE	8-0	60	3/8 CIRCLE	 ADVAPASS	 ADVAPASS	9	12
422	ADVA8732	BLUE	8-0	60	3/8 CIRCLE			6.5	12
423	ADVA8741	BLUE	8-0	60	3/8 CIRCLE	 MICRO VASCULAR	 MICRO VASCULAR	8	12
424	ADVA8741H	BLUE	8-0	60	3/8 CIRCLE	 MICRO VASCULAR	 MICRO VASCULAR	8	36
425	ADVA8753	BLUE	8-0	60	3/8 CIRCLE			9.3	12
426	ADVA1710	BLUE	10-0	30	3/8 CIRCLE			6.5	12
427	ADVA1713	BLUE	10-0	23	STRAIGHT	 MICROPOINT	 MICROPOINT	16	12
428	ADVA8095	BLUE	10-0	75	STRAIGHT			16	12
429	ADVA8155	BLUE	10-0	30	1/2 CIRCLE & STRAIGHT CUT			16	12

 Round Bodied

 Spatulated Micropoint

 Straight Cut

 Taper Cut

 Cutting


Product Description

Composition	Polypropylene
Coating	None
Sizes	U.S.P. Size 2 to 10-0
Colour	Blue
Absorption Profile	Permanent
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied, Diamond Point & Sternum Reverse Cutting

Tensile Strength

86% more knot pull tensile strength Vs. U.S.P. limit ensures permanent wound support.


Knot pull tensile strength (KgF)¹



REFERENCE _____

1 : Data on file with AMS - Bench test done on Polypropylene U.S.P. size 1, 40 mm Round bodied needles

Monofilament Material



The only Polypropylene with an advanced dispensing system

ADVATRAY™ - the confidence to
suture with ease

ADVA**TRAY**™ a unique
dispensing system
designed to deliver ease
of suture dispensing
with minimal memory in
the suture strand⁷



Minimizing Suture Memory...
Maximizing Knot Security

**Enhanced robust packaging to ensure seal strength
integrity and sterility assurance throughout shelf life**

Advancing the Science of Cardiovascular Needles

Premium Needles for consistent Superior Performance

	ADVAPRIME SPECIALLY ENGINEERED DESIGN	Advanced point geometry and intelligent engineering assures steady, superior control and smoothness
	ADVAGUARD BLUNT POINT	ADVAGUARD is engineered to diminish the chance of needle stick injuries, delivering both optimal strength and ductility with increased precision
	ADVABLACK VISIBILITY ENABLING CONFIDENCE	ADVABLACK needles have been designed to offer outstanding visibility. Specially crafted for enhanced assurance and confidence in the operating zone

ADVAPASS advanced coating technology offers consistency in ease of penetration and maintains needle sharpness and performance from start to finish pass after pass

Minimal tissue trauma enables faster recovery



Advanced MedTech Solutions

P-21-22, 25-26 & 34-35, GIDC, Manjusar, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India
T : (91) 2667 671302 | amsltd.com | info@amsltd.com

Singapore Office : 111, Somerset Road, 07-08, 111 Somerset, Singapore - 238164.
T: (65) 636910 80 / 89 | E: info@amsltd.com | amsltd.com