


N&S has been designed to replace, through surgery, the natural testicle in shape and size in patients who had suffered from testicular agenesis, or when the testicle has been removed through surgery due to several pathologies.

N&S testicular prosthesis is a single use sterile permanent implant made of elastomeric of low-hardness silicone, which is very similar to the natural testicle.



N&S is packaged within a sterile double blister and then in a cardboard box.

CONTENTS

There are different models of N&S prosthesis according to the size. They are:

	<i>T-Small</i>		<i>T-Medium</i>		<i>T-Large</i>		<i>T-XLarge</i>	
	Reference	Qty	Reference	Qty	Reference	Qty	Reference	Qty
Testicular prosthesis	T-Small	1 (one)	T-Medium	1 (one)	T-Large	1 (one)	T-XLarge	1 (one)

1. TECHNICAL SPECIFICATIONS

<i>REFERENCE</i> <i>SPECIFICATION</i>	T-Small	T-Medium	T-Large	T-XLarge
Diameter A	31 mm	37 mm	42 mm	47 mm
Diameter B	23 mm	28 mm	32 mm	34 mm
$R_n = (B/2)^2 / (A/2)$	8,5 mm	10,6 mm	12,2 mm	13,5 mm
$R_h = (A/2)^2 / (B/2)$	20,9 mm	24,4 mm	27,6 mm	32,5 mm
Material	Silicone			
Consistence	Solid			
Hardness (EI-02-01)	35 – 65 Shore 00			
Method of Sterilization	Ethylene Oxide			
Shelf Life	3 years			
Contact time	Permanent more than 30 days			

